

Be selective



## SAMPLING & SAMPLE PREPARATION CATALOG

**Solid Phase Extraction  
SPE Disks**

**96 well plates SPE & Spin SPE**

**Microelution & SPE Tips**

**Passive sampling  
and associated products**



Order on  
[www.affinisep.com](http://www.affinisep.com)



## **Dear Customer**

*We are pleased to offer you our new catalog for **SAMPLING & SAMPLE PREPARATION** containing innovative solutions to help you to achieve your goals.*

*AFFINISEP is expert in making consumables for sample preparation and extraction techniques. We save you time by simplifying your workflows with smart application-specific kits, ready to use pre-developed methods, and the highest level of scientific support for your selective extraction needs. Our R&D team is fully committed to developing high quality products that help you to get **accurate** information as **quickly** as possible;*

*We provide knowledge and expertise to produce superior off-the-shelf consumables and custom resins. Our manufacturing organization supported by highly skilled sales and support team is at your disposal to satisfy you.*

*And particularly, we are committed to providing the best technical support possible. Our in-house scientific team will work with you to develop personalized solutions to fit your specific application. For technical inquiries, please contact us by email: **tech.support@affinisep.com***

*We are very thankful to hear all your feedbacks about our products, protocols and customer services by email to: **contact@affinisep.com***

*Your satisfaction is our engine. Thank you to trust us!*

*Kaynoush Naraghi  
CEO*



# TABLE OF CONTENTS

<b>PRODUCT SELECTION GUIDE</b>	<b>8</b>
<b>SOLID PHASE EXTRACTION SOLUTIONS</b>	<b>9</b>
<b>PRODUCTS PRESENTATION</b>	
<b>SOLID PHASE EXTRACTION</b>	
<b>FORMATS</b>	<b>13</b>
<b>SOLID PHASE EXTRACTION SOLUTIONS – Proteomics and biological applications</b>	<b>16</b>
AttractSPE® Disks Technology for microelution	18
Focus on AttractSPE® Disks proteomics tools	20
AttractSPE® Disks chemistries	22
AttractSPE® Tips – Stagetips	23
AttractSPE® Disks 96 plate for microelution	25
AttractSPE® Disks 96 plate – 1mL	26
AttractSPE® Spin column for microextraction	27
AttractSPE® Disks Cartridges	30
<b>SOLID PHASE EXTRACTION SOLUTIONS – For simple to very complex matrices</b>	<b>32</b>
<b>Quick look-up guide</b>	
<b>READY TO USE SPECIFIC KITS</b>	
<b>AFFINIMIP®SPE (MOLECULARLY IMPRINTED POLYMER) AND DEDICATED APPLICATIONS</b>	<b>36</b>
<b>AFFINIMIP® SPE – CONCEPTS</b>	<b>37</b>
<b>AFFINIMIP® SPE FOR ANALYTICAL PURPOSE</b>	<b>39</b>
<b>AFFINIMIP® SPE FOR MYCOTOXINS</b>	<b>40</b>
AFFINIMIP®SPE <b>Multimyc</b> <b>LCMSMS</b> for Multimycotoxins analyses	41
AFFINIMIP®SPE <b>FumoZON</b> for Fumonisin and Zearalenone (simultaneously)	42
AFFINIMIP®SPE <b>Patulin</b>	43
AFFINIMIP®SPE <b>Ochratoxin A</b>	44
AFFINIMIP®SPE <b>Deoxynivalenol</b>	45
AFFINIMIP®SPE <b>Zearalenone</b>	46

# TABLE OF CONTENTS

AFFINIMIP®SPE <b>Glyphosate</b> for glyphosate, AMPA, glufosinate...	47
AFFINIMIP®SPE <b>Picolinic herbicides</b> for Picloram, Aminopyralid, Clopyralid	48
AFFINIMIP®SPE <b>PAHs</b>	49
AFFINIMIP®SPE <b>Tetracyclines</b>	50
AFFINIMIP®SPE <b>Chloramphenicol</b>	51
AFFINIMIP®SPE <b>Aminoglycosides</b>	52
AFFINIMIP®SPE <b>Estrogens</b>	53
AFFINIMIP®SPE <b>Bisphenols</b>	54
AFFINIMIP®SPE <b>Zeranol residues</b>	55
AFFINIMIP®SPE <b>Phenolics</b>	56
AFFINIMIP®SPE <b>NNAL</b>	57
AFFINIMIP®SPE <b>Amphetamines</b>	58
AFFINIMIP®SPE <b>Catecholamines</b>	59
AFFINIMIP®SPE <b>Metanephrines</b>	60
AFFINIMIP®SPE <b>Tamoxifen</b>	61
AFFINIMIP®SPE <b>Cannabis</b>	62
AttractSPE® <b>PFAS</b> for perfluorinated compounds	63
AttractSPE® <b>Non-target screening</b>	64
AttractSPE® <b>Chelating</b> for metal retention	64
SilactSPE <b>DAU</b> (Drug of abuse use) & SilactSPE <b>Organotins</b>	65
SPE for Polycyclic aromatic Hydrocarbons with AFFINIMIP®SPE <b>PAHs</b> , SilactSPE <b>CN/SiOH</b> & AttractSPE® <b>HLB</b>	66
SPE for the interferences removal – AttractSPE® <b>SAX-HCO<sub>3</sub></b> & AttractSPE® <b>PS-H</b>	
SPE for the interferences removal – AttractSPE® <b>PS-Ag</b> & AttractSPE® <b>PS-Ba</b> & SilactPE® <b>HydroxyApatite</b>	68

## SPE FOR GENERAL METHODS


<b>AttractSPE® POLYMERIC-BASED SPE</b>	69
AttractSPE® <b>HLB</b> - Hydrophilic Lipophilic Balance	72
AttractSPE® <b>WAX</b> - Weak Anion Exchange	73
AttractSPE® <b>WCX</b> - Weak Cation Exchange	74
AttractSPE® <b>SAX</b> - Strong Anion Exchange	75
AttractSPE® <b>SCX</b> - Strong Cation Exchange	76
AttractSPE® <b>DVB</b> - Reversed phase	77
<b>Inorganic based SilactSPE CARTRIDGES</b>	78
SilactSPE <b>C18</b> (end capped) & SilactSPE <b>C18NEC</b> (not end capped)	80
SilactSPE <b>C8</b> & SilactSPE <b>Phenyl</b>	81

# TABLE OF CONTENTS

SilactSPE Silica & SilactSPE Cyano	82
SilactSPE Amine (SiWAX) & SilactSPE PSA	83
SilactSPE SiWCX & SilactSPE SiSCX	84
SilactSPE SiSAX & SilactSPE Carbonate	85
SilactSPE Florisil PR, SilactSPE Dry & SilactSPE Na <sub>2</sub> SO <sub>4</sub> / Florisil	86
SilactSPE Alumina-Acidic, Neutral and Basic	87
AttractSPE® Carbon based SPE with AttractSPE® carbon, AttractSPE® carbon/Amine & AttractSPE® carbon/PSA	88
<b>AttractSPE® On-Line SPE</b>	<b>89</b>
<b>Filtration &amp; Passthrough methods</b>	
<b>AttractSPE® SLE, AttractSPE® LipRem, PROTEINS &amp; LIPID REMOVAL; AttractFiltr</b>	
SPE for the removal of proteins or lipids – AttractSPE® LipRem & SilactSPE double fritted or single Fritted	91 92
AttractSPE® SLE for Supported Liquid Extraction	93
AttractFiltr for membrane filter cartridge	94
<b>QUECHERS AND EXTRACTION SALTS</b>	<b>95</b>
<b>WATER MONITORING SPE DISKS &amp; PASSIVE SAMPLING</b>	<b>99</b>
<b>AttractSPE® DISKS for environmental applications</b>	<b>100</b>
AttractSPE® Disks Environment	103
AttractSPE® Disks Manifolds and accessories	104
<b>PASSIVE SAMPLING SOLUTIONS – POCIS, SPATT, Chemcatchers and other passive samplers using disks, silicone rubber strips</b>	<b>105</b>
POCIS product list	109
POCIS for groundwater & Disks passive samplers product list	110
SPATT, Silicone rubber and accessories product list	111
<b>SPE ACCESSORIES</b>	<b>112</b>
AttractSPE® Positive Pressure Manifold	113
Vacuum SPE Manifold and ACCESSORIES	114
<b>EXAMPLES OF SPE APPLICATIONS</b>	<b>115</b>
<b>CUSTOM-MADE PRODUCTS &amp; PRIVATE LABEL SERVICES</b>	<b>122</b>
Custom-made products & Services	123
Project development	124
Private label services	125
<b>ORDERING INFORMATION</b>	<b>126</b>

# PRODUCT SELECTION GUIDE

Sampling and sample preparation are the key steps in trace analysis for analytical chemist. As specialist in this field, AFFINISEP supplies a complete range of solutions based on Solid Phase Extraction (SPE), passive sampling and filtration processes for automated, high-throughput or manual sample preparation.



A multitude of chemical phases and formats are available for various applications. We offer a comprehensive range of SPE to give you all elements to face the increasingly complex and diverse sample preparation challenges by:

- **Simplify data analysis by removing interferences**
- **Increase sensitivity and reliability by enrichment of the analyte**
- **Obtain high and reproducible recovery yields from complex samples**

# Solid Phase Extraction Solutions



2L -  
500mL

- AttractSPE® Disk
- Passive samplers: POCIS, SPATT, Chemcatcher



Water Analysis



200mL  
- 1mL

- AFFINIMIP® SPE
- AttractSPE®
- SilactSPE®
- Qcleanup® Quechers
- AttractFiltration® Filtration
- AttractSPE® SLE



Environmental &  
Food Analysis



1000µL  
- 10µL

- AttractSPE® Disks Tips – Stagetips
- AttractSPE® Disks 96 well-plate
- AttractSPE® Disks Spin SPE
- AttractSPE® Disks SPE cartridges



Protein, Peptides  
& Bioanalysis

Order on  
[www.affinisep.com](http://www.affinisep.com)

# PRODUCTS PRESENTATION

## Solid phase extraction (SPE)

**AFFINIMIP®SPE, AttractSPE® and SilactSPE** are a powerful techniques to provide a rapid and selective sample clean-up, a high recovery and the concentration necessary for accurate quantitative analysis. Usually available as a SPE cartridge format with flange or without flange (Flangeless for some automate applications)

## SPE Disks

**AttractSPE®Disks** are Solid Phase Extraction Disks for the extraction of a broad range of contaminants. **AttractSPE®Disks** are thin and uniform SPE disks for retention of targeted analytes without any breakthrough. **Attract SPE®Disks** have been diversified in two very different fields requiring their own specifications: Environment and Biomolecular analysis.

## Passive sampling for microcontaminants or biotoxins

Passive sampling enables the monitoring of contaminants in water for a short to long period (with an average field deployment of one month) for which no power, maintenance and supervision is required. An average of the concentration of collected contaminants (pesticides, drugs residues, biotoxins, PAHs, Glyphosate & AMPA...) is measured in the laboratory. AFFINISEP provides a wide range of solutions for passive sampling for **POCIS, SPATT, Chemcatcher** for the sampling of contaminants such as pesticides, drugs or biotoxins in water.

## Filtration with Membranes or with frits

**AttractFiltra** is a filtration cartridge or 96-well plate based on the use of membrane for sample particles filtration using vacuum manifold, gravity or SPE automates. A broad range of membranes is available and can enable a broad range of sample filtration. Double fritted cartridges are also available for protein precipitation.

## Dispersive SPE - Quechers

**Qcleanup** dispersive SPE products are a mixture of Magnesium Sulfate, PSA, Black carbon or C18 used by QuEChERS method for the clean-up of fruits and vegetables during multiresidues pesticides analyses.

## Supported Liquid Extraction (SLE)

**AttractSPE®SLE** is Supported Liquid Extraction (*a.k.a* SLE) , an alternative to LLE to pass from an aqueous media to an organic media without emulsion formation. SLE is used to transfer hydrophobic molecules from water to an organic solvent. In addition, this process can be used for the removal of proteins or lipids.

## Microextractions

For proteomics or genomics analysis, sample size is often a few microliters. However, sample preparation is required and sample preparation tools must work for such small sample size.

Some sample preparation tools are presented below. They are based on the use of **AttractSPE® Disks** as sorbent. Thanks to these densely packed beads embedded in the disk, the volume of elution is reduced in comparison to classical packed powder.

### **AttractSPE® Tips – Stagetips**

μSPE column designed by immobilizing a uniform disk inside a pipette tips (Stage-tips).

### **AttractSPE® Spin SPE**

SPE column designed by immobilizing a uniform disk inside a microcentrifuge SPE tube.

### **AttractSPE® Disks 96 well-plate SPE**

96 SPE well plate designed by immobilizing a uniform disk.

### **AttractSPE® Disks SPE cartridges**

SPE cartridges designed by immobilizing a uniform disk.

# SOLID PHASE EXTRACTION

## SPE PROCEDURE STEPS


**1- Sample preparation:** This step is necessary to obtain a loading solution compatible with retention condition in the SPE column. In case of solid matrices, this solution must also extract the compounds of interest from these matrices.

**2- Conditioning:** The SPE cartridges are conditioned with the appropriated solvents to fully soak sorbent and enable further interactions between the analytes and the sorbent.

**3- Loading:** the percolation solution is loaded through the SPE cartridge. The analyte must be retained in the column as well as unwanted compounds.

**4- Washing:** interferences and unwanted compounds are washed off by using appropriate solvents.

**5- Elution:** The desired analyte is extracted from the SPE cartridge.



### Normal phase

- Based on polar-polar or dipole-dipole interactions between the analyte and a non organic phase like silica.

### Reversed phase

- Based on non polar- non polar interactions and Van der Waals dispersive forces. The sorbent is hydrophobic like polymeric sorbent modified silica-based sorbent.

### Ion-exchange

- Uses electrostatic interaction between a charged sorbent and the ionic analyte. The sorbent is charged with the opposite charge of the analyte.

### Mixed-mode sorbents

- Interact through reversed phase and ion exchange retention mechanisms. Available as a polymeric sorbent (AttractSPE® SAX, WAX, WCX or SCX) or as SilactSPE .

### Imprinted

- Highly selective based on forme and interaction of one molecule or a family of molecules.

# Formats

The choice of format is essential for a successful analysis.

Our products are available under different formats. We can provide as well other formats under request.

If you need some support to identify the right format, please contact our technical support.

For any question, please contact us:

[tech.support@affinisep.com](mailto:tech.support@affinisep.com)

## Open Cartridge



Formats: 1mL, 3mL, 6mL,  
15mL, 20mL, 60mL  
With flange or tabless  
Materials: Polypropylene,  
glass (6mL)  
Frits: Polyethylene, PTFE ,  
Glass fiber  
Sorbents: powder or disk  
*luer compatible*

## Reversible Cartridge

Formats: 0,7mL, 2mL  
Material: Polypropylene  
Frit: Polyethylene  
*Luer compatible*



## On-line SPE Cartridge

I.D: 2,1 and 4,6mm  
Length: 20mm



## LRC Cartridge



*Luer compatible*  
Formats: 10mL

Material: Polypropylene  
Frit: Polyethylene

## Cartridge for Automates

Formats: 1mL, 3mL, 6mL  
Material: Polypropylene  
Frit: Polyethylene





1000 $\mu$ l - 10 $\mu$ l

- AttractSPE<sup>®</sup> Tips – Stagetips
- AttractSPE<sup>®</sup> Spin SPE
- AttractSPE<sup>®</sup> Disks 96 well-plate SPE
- AttractSPE<sup>®</sup> Disks SPE cartridges

## Proteomics and biological applications

# AttractSPE® Disks

Proteomics, Biomarker discoveries  
and Biological applications

Spinnable, Automatable  
High throughput HTS  
Micro-SPE for MicroElution

Proteins, Peptides, DNA, Small molecules, ...  
Purification,  
Desalting  
Fractionation  
Bioanalysis



Order on  
[www.affinisep.com](http://www.affinisep.com)

## AttractSPE® Disks Technology for microelution

AttractSPE®Disks are **thin, dense, soft and uniform membranes** based chromatography for extraction/separation, purification and concentration of analyte molecules.

**Thanks to their unique advantage,** AttractSPE®Disks are useful for purification of Very Small Sample Volumes in Proteomics, Genomics, Metabolomics, Bioanalysis, Biomarker discoveries and Biological applications. They are applied for Spinnable, Automatable, High throughput microelution and nano elution.

AttractSPE®Disks offer outstanding sample preparation efficiency and reproducibility of results. Since the diffusion distance between particles is minimized, adsorption is more efficient, and extraction can be accomplished using **very low sorbent mass**.

These properties are giving to AttractSPE®Disks a significant improvement of mass transfer kinetics compared to traditional packed SPE particles. As a monolith disk, AttractSPE®Disks are self stand and require no frits for immobilizing the column bed (unlike traditional SPE products) allowing 100% recovery of the original sample volume.

### AttractSPE®Disks Advantages

- No need to frits or filters
- Reduced dead volume
- Small elution volumes
- High sample recovery
- Reduced time for eluate evaporation
- Higher throughput
- Channeling effects eliminated
- Excellent reproducibility
- Concentration of the sample



# AttractSPE® Disks Technology for microelution

AttractSPE® Disks can be used to miniaturize SPE for small operating volume of fluid samples. **4 formats** are **available** for microelution. Thanks to the use of AttractSPE® Disks **for all formats**, then **the change of format or the scale up of the process is easy**.

Capacity of molecules of interest increases

## □ AttractSPE® Tips – Stagetips and AttractSPE® Disks 96 well-Plate for microelution

μSPE columns designed by immobilizing a uniform disk inside a pipette tips (Stagetips per unit or as 96 well plate)



## □ AttractSPE® Disks 96 well-Plate

96 SPE well plate designed by immobilizing a uniform disk



## □ AttractSPE® Spin Columns

SPE column designed by immobilizing a uniform disk inside a centrifuge SPE tube



## □ AttractSPE® Disks Cartridges

SPE cartridges designed by immobilizing a uniform disk



Concentration of molecules of interest increases

## Focus on AttractSPE® Disks proteomic tools

AttractSPE® **Tips**, **Spin**, **Disks 96 plate** and **Disks cartridges** are tools for proteomic applications (fractionation and desalting of peptides/proteins) and for bioanalysis (small molecules cleanup). In order to give the most exhaustive applications, they are proposed with a large variety of sorbents (HLB, RPS, SDB, SCX, SAX, C8, C18...).



AttractSPE®  
**Tips**



AttractSPE®  
**Spin**



AttractSPE® **Disks 96**  
**plate for microelution**




AttractSPE®  
**Disks cartridges**

### Advantages of AttractSPE® Disks tools

- Removes interfering contaminants—significantly reduces signal suppression and improves signal-to-noise ratios and sequence coverage
- Simplifies optimization—processing yields high-quality spectra and is effective for a variety of reverse-phase-compatible contaminants
- Robust—works with a wide variety of load volumes and concentrations; no need to reduce sample volume before application
- Convenient—easy to handle and requires no special equipment to process multiple samples simultaneously (unlike tip-driven systems that require one sample to be processed at a time)
- Sensitive—special C18 resins allows excellent recovery percentages, even at low (sub-picomole) sample loads

### Capacity

Each product is available with different capacities, which are designated T1, T2 or T3 in the product designations. This designation refers to the layer thickness of SPE tips and the capacity. For a more reproducible product, each tool contains only **ONE** layer. **References below (p24 to p30) will have **XX** indication which has to be replaced by the capacity you require (T1, T2 or T3) on the reference.**



*Capacity for 200µL tips*

- T1 : capacity up to 15µg
- T2 : capacity up to 30µg
- T3 : capacity up to 50µg

### AttractSPE® Disks sorbents

- A broad variety of sorbents for each required applications
- Various formats: disks, spins, 96 SPE plates, cartridges
- One sorbent - several layers for increased capacity
- Several sorbents - stacking for complex applications
- Disks used as filter for application requiring beads

# AttractSPE® Disks chemistries



Sorbents for SPE Disks for biomolecular applications	Compatible with analytical methods
<ul style="list-style-type: none"> <li>- <b>C18 end-capped (C18 EC)</b> for the most hydrophobic compounds/peptides</li> <li>- <b>C18 Not end-capped (C18 NEC)</b> with the best combination - large spectrum of interaction - capacity</li> <li>- <b>C18</b> for the wider interaction spectrum (hydrophobic - hydrophilic) with compounds / peptides</li> </ul>	<ul style="list-style-type: none"> <li>- Desalting of peptides; fractionation of peptides at acidic and neutral pH</li> <li>- Drug extraction in biological samples,</li> </ul>
<b>C8</b>	Desalting of large peptides and proteins; Usage as frits to retain beads in a tip
<b>Silica</b>	Purification of DNA
<b>C4</b>	Desalting of large peptides and proteins
<b>SDB</b> <i>a.k.a</i> PS-DVB	Fractionation of peptides at basic pH
<b>HLB:</b> SDB with hydrophilic moieties	Fractionation of peptides Extraction of small molecules (drugs) in biological fluids
<b>RPS:</b> Sulfonic modified SDB sorbent	Desalting of peptides; fractionation of peptides
<b>SAX</b> : Anion exchange SDB	Fractionation of peptides by salt or pH steps
<b>SCX</b> : Cation exchange SDB	Fractionation of peptides by salt or pH steps

## AttractSPE® Tips – StageTips

**AttractSPE® Tips are spinnable and automatable Tips for high throughput useful for peptide desalting, proteomics, small molecules screening and biomarker discoveries and biological applications**

Tip Volume (µL)	Resin Volume (µL)	Resin:	Application
<ul style="list-style-type: none"><li>•200</li><li>•300</li><li>•500</li><li>•1000</li></ul>	<ul style="list-style-type: none"><li>•5</li><li>•10</li><li>•12</li><li>•20</li><li>•40</li><li>•80</li><li>•160</li><li>•200</li></ul>	<ul style="list-style-type: none"><li>•C18 300Å</li><li>•C18 100Å</li><li>•C8 100Å</li><li>•C4 300Å</li><li>•Strong Anion</li><li>•Strong Cation</li><li>•Weak Anion</li><li>•Weak Cation</li><li>•HLB</li><li>•RPS</li></ul>	<ul style="list-style-type: none"><li>•Desalting</li><li>•PCR Purification kit</li><li>•Genomic Purification</li><li>•Fractionation</li><li>•Plasmid Purification</li><li>•Bioanalysis...</li></ul>

To prepare high quality peptide samples for LC-MS, it is very important to ensure the overall quality of shotgun proteomics experiments. Peptide samples collected after digestion usually need to be cleaned to remove salts, possible gel pieces (for in-gel digested samples) or particles (for in-solution digested samples), which otherwise will damage the LC switching valves or clog the columns.



# AttractSPE® Tips – Stage Tips

## Advantages

- Load your sample on AttractSPE® Disks Tips for desalting or fractionate peptides/ proteins – bioanalysis of small molecules
- Several sorbents based Stage-tips and stacking
- Available as 10, 20, 100, 200µL, 1mL




Designation	Description	Reference – 200µl - 96/pk **	Reference –1mL - 96/pk**
<b>AttractSPE® Tips C18</b>	C18 membrane, 96/pk	Tips-C18.XX.200.96	Tips-C18.XX.1000.96
<b>AttractSPE® Tips C18 EC</b>	End-capped C18 membrane, 96/pk	Tips-C18EC.XX.200.96	Tips-C18EC.XX.1000.96
<b>AttractSPE® Tips C18 NEC</b>	Not end-capped C18 membrane, 96/pk	Tips-C18NEC.XX.200.96	Tips-C18NEC.XX.1000.96
<b>AttractSPE® Tips C8</b>	C8 membrane, 96/pk	Tips-C8.XX.200.96	Tips-C8.XX.1000.96
<b>AttractSPE® Tips SDB</b>	PS-DVB membrane, 96/pk	Tips-DVB.XX.200.96	Tips-DVB.XX.1000.96
<b>AttractSPE® Tips RPS</b>	<b>Modified DVB</b> membrane, 96/pk	Tips-RPS-M.XX.200.96	Tips-RPS-M.XX.1000.96
<b>AttractSPE® Tips HLB</b>	HLB membrane, 96/pk	Tips-HLB.XX.200.96	Tips-HLB.XX.1000.96
<b>AttractSPE® Tips SAX</b>	SAX membrane, 96/pk	Tips-SAX.XX.200.96	Tips-SAX.XX.1000.96
<b>AttractSPE® Tips SCX</b>	SCX membrane, 96/pk	Tips-SCX.XX.200.96	Tips-SCX.XX.1000.96
<b>AttractSPE® Tips C18-SCX</b>	Stacking C18 & SCX membranes, 96/pk	Tips-C18-SCX.XX.200.96	Tips-C18-SCX.XX.1000.96
<b>AttractSPE® Tips C18-SCX-C18</b>	Stacking C18&SCX&C18 membranes, 96/pk	Tips-C18-SCX-C18.XX.200.96	Tips-C18-SCX-C18.XX.1000.96
<b>AttractSPE® Tips Silica</b>	Silica membranes, 96/pk	Tips-Si.XX.200.96	Tips-Si.XX.1000.96

\*\* Replace XX by the capacity T1, T2 or T3 depending on your needs.  
e.g.: Tips-C18.T3.200.96

## AttractSPE® Disks 96 plate for microelution

AttractSPE® Disks 96 Plate for microelution is a 96 plate with wells similar to stagetips for a high throughput clean-up requiring microelution volumes.



Designation	Description	Reference – 1/pk **
AttractSPE® Disks 96 plate <b>C18</b>	C18 membrane, 1/pk	μ96W-C18.XX.1
AttractSPE® Disks 96 plate <b>C18 EC</b>	End-capped C18 membrane, 1/pk	μ96W-C18EC.XX.1
AttractSPE® Disks 96 plate <b>C18 NEC</b>	Not end-capped C18 membrane, 1/pk	μ96W-C18NEC.XX.1
AttractSPE® Disks 96 plate <b>C8</b>	C8 membrane, 1/pk	μ96W-C8.XX.1
AttractSPE® Disks 96 plate <b>SDB</b>	PS-DVB membrane, 1/pk	μ96W-DVB.XX.1
AttractSPE® Disks 96 plate <b>RPS</b>	<b>Sulfonated Modified DVB</b> membrane, 1/pk	μ96W-RPS-M.XX.1
AttractSPE® Disks 96 plate <b>HLB</b>	HLB membrane, 1/pk	μ96W-HLB.XX.1
AttractSPE® Disks 96 plate <b>SAX</b>	SAX membrane, 1/pk	μ96W-SAX.XX.1
AttractSPE® Disks 96 plate <b>SCX</b>	SCX membrane, 1/pk	μ96W-SCX.XX.1
AttractSPE® Disks 96 plate <b>WCX</b>	WCX membrane, 1/pk	μ96W-WCX.XX.1
AttractSPE® Disks 96 plate <b>C18-SCX</b>	Stacking C18 & SCX membranes, 1/pk	μ96W-C18-SCX.XX.1
AttractSPE® Disks 96 plate <b>C18-SCX-C18</b>	Stacking C18 & SCX & C18 membranes, 1/pk	μ96W-C18-SCX-C18.XX.1
AttractSPE® Disks 96 plate <b>Silica</b>	Silica membranes, 1/pk	μ96W-Si.XX.1

\*\*Replace XX by the capacity (T1, T2, T3) depending on your needs.

## AttractSPE® Disks 96 plate – 1mL

**AttractSPE® Disks 96 Plate** is a plate with 96 SPE wells of 1mL containing immobilized SPE disks that enables a high throughput clean-up with the simultaneous preparation of 96 samples. Thanks to this small sorbent amount and a high efficiency, almost 100% of the original sample is recovered.

### Advantages

- All sorbents available on catalog or on demand
- Several sorbent weights available.
- Easy handling with automates or liquid handling robots - spinnable
- AttractSPE® Manifold for 96 wellPlate format or AttractSPE® Positive pressure Manifold for 96 well Plate




Designation	Description	Reference – 1/pk**
<b>AttractSPE® Disks 96 plate C18</b>	C18 membrane, 1/pk	96W-C18.XX.1
<b>AttractSPE® Disks 96 plate C18 EC</b>	End-capped C18 membrane, 1/pk	96W-C18EC.XX.1
<b>AttractSPE® Disks 96 plate C18 NEC</b>	Not end-capped C18 membrane, 1/pk	96W-C18NEC.XX.1
<b>AttractSPE® Disks 96 plate C8</b>	C8 membrane, 1/pk	96W-C8.XX.1
<b>AttractSPE® Disks 96 plate SDB</b>	PS-DVB membrane, 1/pk	96W-DVB.XX.1
<b>AttractSPE® Disks 96 plate RPS</b>	<b>Modified</b> DVB membrane, 1/pk	96W-RPS-M.XX.1
<b>AttractSPE® Disks 96 plate HLB</b>	HLB membrane, 1/pk	96W-HLB.XX.1
<b>AttractSPE® Disks 96 plate SAX</b>	SAX membrane, 1/pk	96W-SAX.XX.1
<b>AttractSPE® Disks 96 plate SCX</b>	SCX membrane, 1/pk	96W-SCX.XX.1
<b>AttractSPE® Disks 96 plate C18-SCX</b>	Stacking C18 & SCX membranes, 1/pk	96W-C18-SCX.XX.1
<b>AttractSPE® Disks 96 plate C18-SCX-C18</b>	Stacking C18 & SCX & C18 membranes, 1/pk	96W-C18-SCX-C18.XX.1
<b>AttractSPE® Disks 96 plate Silica</b>	Silica membranes, 1/pk	96W-Si.XX.1

\*\*Replace XX by the capacity (T1, T2, T3) depending on your needs. e.g.: 96W-C18. T1.1

## AttractSPE® Spin Column for microelution

**AttractSPE® Spin Column** is an SPE column created by immobilizing a monolithic disk inside a microcentrifuge SPE tube. Thanks to its self-stand, the monolith disk requires no frits for immobilizing the column bed (unlike traditional SPE products), which allows essentially 100% recovery of the original sample volume after a couple of brief centrifugations.



# AttractSPE® Spin Column

## Advantages

- 4 spin size formats: micro and minispin, 15mL and 50mL spin tubes
- Fast and easy extraction process by centrifugation
- High throughput purification
- Easy scale up of the process



Mini spin

Micro spin

Designation	Description	Reference micro spin – 96/pk**	Reference mini spin – 96/pk**
<b>AttractSPE® Spin C18</b>	C18 membrane, 96/pk	μSpin-C18.XX.96	Spin-C18.XX.96
<b>AttractSPE® Spin C18 EC</b>	End-capped C18 membrane, 96/pk	μSpin-C18EC.XX.96	Spin-C18EC.XX.96
<b>AttractSPE® Spin C18 NEC</b>	Not end-capped C18 membrane, 96/pk	μSpin-C18NEC.XX.96	Spin-C18NEC.XX.96
<b>AttractSPE® Spin C8</b>	C8 membrane, 96/pk	μSpin-C8.XX.96	Spin-C8.XX.96
<b>AttractSPE® Spin SDB</b>	PS-DVB membrane, 96/pk	μSpin-DVB.XX.96	Spin-DVB.XX.96
<b>AttractSPE® Spin RPS</b>	Modified DVB membrane, 96/pk	μSpin-RPS-M.XX.96	Spin-RPS-M.XX.96
<b>AttractSPE® Spin HLB</b>	HLB membrane, 96/pk	μSpin-HLB.XX.96	Spin-HLB.XX.96
<b>AttractSPE® Spin SAX</b>	SAX membrane, 96/pk	μSpin-SAX.XX.96	Spin-SAX.XX.96
<b>AttractSPE® Spin SCX</b>	SCX membrane, 96/pk	μSpin-SCX.XX.96	Spin-SCX.XX.96
<b>AttractSPE® Spin C18-SCX</b>	Stacking C18 & SCX membranes, 96/pk	μSpin-C18-SCX.XX.96	Spin-C18-SCX.XX.96
<b>AttractSPE® Spin C18-SCX-C18</b>	Stacking C18 & SCX & C18 membranes, 96/pk	μSpin-C18-SCX-C18.XX.96	Spin-C18-SCX-C18.XX.96
<b>AttractSPE® Spin Silica</b>	Silica membranes, 96/pk	μSpin-Si.XX.96	Spin-Si.XX.96
<b>Reservoirs for AttractSPE® Spin</b>	25mL- 96/pk	Spin-Res.96	Spin-Res.96

\*\*Replace XX by the capacity (T1, T2, T3) depending on your needs. e.g.: μSpin-C18.T2.96

## AttractSPE® Spin Column

AttractSPE® Spin columns of 15mL and 50mL are the perfect tools to scale up fractionation and desalting of peptides thanks to a high area of SPE disks.



Designation	Description	Reference 15mL spin tube – 50/pk Loading <4mL**	Reference 50mL spin tube – 50/pk Loading <22mL**
<b>AttractSPE® Spin C18</b>	C18 membrane, 50/pk	Spin15-C18.XX.50	Spin50-C18.XX.50
<b>AttractSPE® Spin C18 EC</b>	End-capped C18 membrane, 50/pk	Spin15-C18EC.XX.50	Spin50-C18EC.XX.50
<b>AttractSPE® Spin C18 NEC</b>	Not end-capped C18 membrane, 50/pk	Spin15-C18NEC.XX.50	Spin50-C18NEC.XX.50
<b>AttractSPE® Spin C8</b>	C8 membrane, 50/pk	Spin15-C8.XX.50	Spin50-C8.XX.50
<b>AttractSPE® Spin SDB</b>	PS-DVB membrane, 50/pk	Spin15-DVB.XX.50	Spin50-DVB.XX.50
<b>AttractSPE® Spin RPS</b>	Modified DVB membrane, 50/pk	Spin15-RPS-M.XX.50	Spin50-RPS-M.XX.50
<b>AttractSPE® Spin HLB</b>	HLB membrane, 50/pk	Spin15-HLB.XX.50	Spin50-HLB.XX.50
<b>AttractSPE® Spin SAX</b>	SAX membrane, 50/pk	Spin15-SAX.XX.50	Spin50-SAX.XX.50
<b>AttractSPE® Spin SCX</b>	SCX membrane, 50/pk	Spin15-SCX.XX.50	Spin50-SCX.T1.50
<b>AttractSPE® Spin C18-SCX</b>	Stacking C18 & SCX membranes, 50/pk	Spin15-C18-SCX.XX.50	Spin50-C18-SCX.XX.50
<b>AttractSPE® Spin C18-SCX-C18</b>	Stacking C18 & SCX & C18 membranes, 50/pk	Spin15-C18-SCX-C18.XX.50	Spin50-C18-SCX-C18.XX.50
<b>AttractSPE® Spin Silica</b>	Silica membranes, 50/pk	Spin15-Si.XX.50	Spin50-Si.XX.50

\*\*Replace XX by the capacity (T1, T2, T3) depending on your needs. e.g.: Spin15-C18. T3.50

# AttractSPE® Disks Cartridges

## Advantages

- 3 and 6mL format
- Larger loading volume with a minimal elution volume
- High extraction capacity
- A broad range of sorbents or sorbent combination



Designation	Reference – 1mL – 50/pk**	Reference – 3mL – 50/pk**	Reference – 6mL – 50/pk**
AttractSPE® Disks Cartridge <b>C18</b>	CAR1-C18.XX.50	CAR3-C18.XX.50	CAR6-C18.XX.50
AttractSPE® Disks Cartridge <b>C18 EC</b>	CAR1-C18EC. XX.50	CAR3-C18EC.XX.50	CAR6-C18EC.XX.50
AttractSPE® Disks Cartridge <b>C18 NEC</b>	CAR1-C18NEC. XX.50	CAR3-C18NEC. XX.50	CAR6-C18NEC. XX.50
AttractSPE® Disks Cartridge <b>C8</b>	CAR1-C8.XX.50	CAR3-C8.XX.50	CAR6-C8.XX.50
AttractSPE® Disks Cartridge <b>SDB</b>	CAR1-DVB.XX.50	CAR3-DVB.XX.50	CAR6-DVB.XX.50
AttractSPE® Disks Cartridge <b>RPS</b>	CAR1-RPS- M.XX.50	CAR3-RPS-M.XX.50	CAR6-RPS-M.XX.50
AttractSPE® Disks Cartridge <b>HLB</b>	CAR1-HLB .XX.50	CAR3-HLB .XX.50	CAR6-HLB .XX.50
AttractSPE® Disks Cartridge <b>SAX</b>	CAR1-SAX .XX.50	CAR3-SAX .XX.50	CAR6-SAX .XX.50
AttractSPE® Disks Cartridge <b>SCX</b>	CAR1-SCX.XX.50	CAR3-SCX.T1.50	CAR6-SCX.XX.50
AttractSPE® Disks Cartridge <b>C18-SCX</b>	CAR1-C18- SCX.XX.50	CAR3-C18-SC.XX.50	CAR6-C18- SCX.XX.50
AttractSPE® Disks Cartridge <b>C18-SCX-C18</b>	CAR1-C18-SCX- C18.XX.50	CAR3-C18-SCX- C18.XX.50	CAR6-C18-SCX-C18. XX.50
AttractSPE® Disks Cartridges <b>SDB-SAX</b>	CAR1-DVB- SAX.XX.50	CAR3-DVB- SA.XX.50	CAR6-DVB- SAX.XX.50
AttractSPE® Disks 96 Cartridges <b>Silica</b>	CAR1-Si.XX.50	CAR3-Si.XX.50	CAR6-Si.XX.50

\*\*Replace XX by the capacity (T1, T2, T3) depending on your needs. e.g.: CAR1-C18. T2.50


# Solid Phase Extraction Solutions

200mL – 1mL

- AFFINIMIP® SPE
- AttractSPE®
- SilactSPE
- Qcleanup Quechers
- AttractFiltra Filtration
- AttractSPE® SLE

**For simple to very  
complex matrices**

**AFFINISEP offers a complete range of sorbents for solid phase extraction from very specific to crude clean-up phases, from silica to polymers, from conventional to more sophisticated sorbents for various applications.**



# SOLID PHASE EXTRACTION

AFFINISEP offers a complete range of chemical phases for solid phase extraction from very specific to crude clean-up phases, from silica to polymers, from conventional to more sophisticated one for various applications.

## **AFFINIMIP®SPE**

**AFFINIMIP®SPE** is a selective solid phase extraction based on Molecularly Imprinted Polymers (MIP) that combines the advantages of immune-affinity columns regarding the selectivity and those of polymeric SPE sorbents in terms of costs and robustness. Developed for complex extraction, an instruction sheet gives all required information to successfully carry out the analysis, including the protocol. No method development is required. The protocol is easy, simple and fast. The chromatographic analysis is fast thanks to the affinity between the analyte and **AFFINIMIP®SPE**. Methods are developed using the most common matrix containing the analyte. These matrices can be as diverse as aqueous matrices (food, feed and biological matrices), fully organic matrices (oil), powdery (milk) or cereals grains. If your sample preparation is complex or gives unreliable results, the development of an **AFFINIMIP® SPE** can be a solution.

## **AttractSPE® Polymeric sorbents**

**AttractSPE®** are the last generation of polymeric SPE sorbent. They are crosslinked polymer particles bearing various chemical functionalities.

These chemical phases provide all the advantages of polymeric material, including robustness (to solvent, pH and T), simplified method development, wide applicability and not affected by drying out. The combination of the water-wettable optimised surface chemistry, high surface area and pH stability ensures high reproducible recoveries for a wide range of analytes. Thanks to their capacity, for the same application, less sorbent and less solvent are required than for silica-based sorbent and so, their added value is particularly obvious for trace analysis and for small volume sample.

## **SilactSPE inorganic-based chemical phases**

The first chemical phases historically used for SPE applications, Silica, Alumine-based sorbents offer a broad range of chemically modified sorbents. This chemistry goes from very polar sorbent to hydrophobic sorbent (end-capped saturated hydrocarbon) passing through intermediate polarity (for instance, amino). These sorbents are sensitive to extreme pH. They are adequate for non-polar analytes in simple matrices.

## QUICK LOOK-UP GUIDE


The choice of the chemical phase is crucial for a successful sample pretreatment. When an **AFFINIMIP® SPE** exists for the analyte, the choice is very simple and the problem is solved. The protocol is simple and ready to use. **AFFINIMIP® SPE** have been developed for various trace analyses in complex matrices and for problematic analysis with common SPE sorbents.

When no **AFFINIMIP® SPE** product is available, you can check the **application notebook** available on our website to use the appropriate sample pretreatment.

In other cases, the sorbent chemistry must be chosen very carefully and each step of the protocol must be optimized. This choice is based on analyte and matrix nature, on the loading condition and on the planned retention mechanism with the sorbent.

The following scheme gives a brief guideline to select the family of sorbents.

For specific applications such as trace analysis or complex matrices, conventional SPE sorbents may give very complex analysis or unsatisfying results. The development of an **AFFINIMIP® SPE** sorbent for this application is a solution. We provide as well services for method development.



# AFFINIMIP® SPE and dedicated applications

## Selective Solid Phase Extraction

Molecularly Imprinted Polymers for the Selective  
Extraction of Trace Analytes from Complex Matrices






## Selective Solid Phase Extraction


### Molecularly Imprinted Polymers for the Selective Extraction of Trace Analytes from Complex Matrices

#### New Extraction Phase Based on Molecularly Imprinted Polymers (MIPs)

MIPs are polymers with shape «memory» and functional groups affine to a template molecule. Using an imprinting process, AFFINISEP designs these materials in order to recognize selectively a target molecule, even in the presence of compounds with structure and functionality similar to the template.



#### Advantages of SPE based molecularly imprinted polymers



# AFFINIMIP® SPE FOR ANALYTICAL PURPOSE

**AFFINIMIP® SPE** is a selective solid phase extraction based on Molecularly Imprinted Polymers (MIP). It combines the advantages of immune-affinity columns regarding the selectivity and of a classic Solid Phase Extraction (SPE) in terms of robustness and costs.




## AFFINIMIP® SPE Selective Extraction Cartridges

### Perfect clean-up system for trace analysis

Thanks to the selectivity of **AFFINIMIP®SPE**, stringent washing steps can be applied in order to remove all interferences and thus minimize matrix effects. It also **reduces ion-suppression effects**.

### Minimal or no method development required

A protocol based on three steps (loading, washing and elution) is supplied with **AFFINIMIP®SPE** kits for tested matrices. No extra-equipment than the usual required for SPE experiments is necessary.



**AFFINIMIP® SPE** protocols are as well defined by 3 steps of loading, washing and elution. All steps have been already developed in detail by AFFINISEP and an instruction sheet is supplied with the product.

# AFFINIMIP® SPE FOR ANALYTICAL PURPOSE




## Benefits

- High affinity and selectivity
- Lowered quantification limits
- High reproducibility and repeatability
- Robust materials

The advantages of **AFFINIMIP® SPE** are essential in trace analysis from a complex matrix in food safety, environment, cosmetics, clinical chemistry, pharmaceutical analysis and others.

The SPE protocol is supplied in an instruction sheet for various complex matrices.

For other matrices, please contact our technical support to help you with your application.



### 7 CLEAN-UP PROCEDURE OF ESTROGENS FROM PLASMA:

**7.1 Preparation of the loading solution**

Dilute your plasma sample by 5 with water. For example, in a 1 mL-volumetric flask: add 0.2 mL of plasma and completed with ultrapure water

**7.2 Protocol for the clean-up of Estrogens from plasma**

Step (Flow rate)	AFFINIMIP® SPE Estrogens (30mg/1mL)
Equilibration with (2 drops/s)	<ul style="list-style-type: none"><li>• 1 mL Acetonitrile</li><li>• 1 mL ultrapure water</li><li>• Do not allow the cartridge to dry during conditioning</li></ul>
Loading (L) (1 drop every 2 seconds)	<ul style="list-style-type: none"><li>• 250µL to 1 mL of the loading solution</li></ul>
Washing of interferents (1 drop/s)	<ul style="list-style-type: none"><li>• 1 mL ultrapure water</li><li>• 1 mL of (60/40) ultrapure water/ Acetonitrile (v/v)</li></ul>
Drying	<ul style="list-style-type: none"><li>• Apply vacuum or nitrogen flow through cartridge during <b>30 seconds</b></li></ul>
Elution (E) (1 drop/s)	<ul style="list-style-type: none"><li>• 1 mL Methanol</li></ul>

The elution fraction (E) is evaporated until dryness under nitrogen with a mini-vap evaporator at room temperature (or a centrifugal evaporator). The residue is dissolved in 0.5 mL of mobile phase for further analysis. Alternatively, the elution may be diluted to a known volume by addition of water for further analysis.

Example of an instruction sheet supplied with **AFFINIMIP® SPE**

# MYCOTOXINS ANALYSES

Mycotoxins are toxic secondary metabolites produced by different fungi present in agricultural commodities. They are regulated in food and feed due to nephrotoxic, neurotoxic, carcinogenic, estrogenic, and immunosuppressive effects.

AFFINISEP has developed two sets of products for mycotoxins analyses:

**Multimycotoxins extraction:** Designed for the simultaneous extraction of several mycotoxins which are present in the same matrix prior to LC-MS/MS analyses. These mycotoxins are all present in the same matrix to be analyzed. Their extraction is done all at once by SPE.

AFFINIMIP® SPE **Multimyco LCMSMS** for the analyses of

- Fumonisin
- Aflatoxin
- Ochratoxin A
- T-2 and HT-2
- Zearalenone
- Deoxynivalenol

AFFINIMIP® SPE **FumoZON** for the analyses of

- Fumonisin
- Zearalenone

**Single mycotoxin extraction:** Designed for the analysis of one specific family of mycotoxin:

AFFINIMIP® SPE **Patulin**

AFFINIMIP® SPE **Ochratoxin A**

AFFINIMIP® SPE **Zearalenone**

AFFINIMIP® SPE **Fumonisin**

AFFINIMIP® SPE **Deoxynivalenol**

AFFINIMIP® SPE **Aflatoxine**



## Analytes

- Fumonisin B1/B2
- Aflatoxins B1/B2/G1/G2
- Ochratoxin A
- T-2 and HT-2
- Zearalenone
- Deoxynivalenol

## Matrices

- Wheat, Maize, Cereals, sunflower seeds, ...

## Advantages

- **SIMULTANEOUSLY** clean-up and concentrate the main regulated mycotoxins

## Regulation

In Europe, these mycotoxins are all regulated (see set values for each individual mycotoxin)  
 Aflatoxins - 1881/2006/EC :  
 foodstuffs : 0,1µg/kg Afla B1 for processed cereals for baby food  
 4,0 µg/kg for total aflas in processed cereals

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
AFFINIMIP® SPE Multimyco LCMSMS	3mL	FS118-03-NG
	6mL	FS118-03B-200-NG

*Visit our website to know more about this product and its formats.  
 Custom designed formats are available on demand.*



## Analytes

- Fumonisin B1+B2
- Zearalenone

## Matrices

- Wheat, Maize, Cereal-based baby food,

## Advantages

- **SIMULTANEOUS** analysis by LC-MS detection
- Very simple protocol

## Regulation

**Europe** - 1126/2007/EC: 20µg/kg and 200µg/kg on maize-based babyfood for respectively Zearalenone and Fumonisin

**Codex Alimentarius** - (193-1995): 2000µg Fumonisin B1+B2/kg for whole commodity

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
AFFINIMIP® SPE FumoZON	3mL	FS109-03

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*



## Analytes

- Patulin

## Matrices

- Apple juice (clear & cloudy), Apple and Multifruit puree, Baby food, Cider, Alcohol, Pommeau, Manzella, Dried apple, Blueberry, Tomato Ketchup

## Advantages

- Applicable to several apple derived matrices
- **Unique** extraction method available on the market
- UV or LC-MS detection

## Regulation

### Regulations for apple puree:

**Europe** - 1881/2006/EC : 25µg/kg; 10µg/kg for infants and young children

**USA (FDA CPG Sec.510.150)** : 50µg/kg

### Regulations for apple juice :

**Europe** - 1881/2006/EC : 50µg/kg; 10µg/kg for infants and young children

**Codex alimentarius:** (193-1995): 50µg/kg apple juice

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE Patulin</b>	3mL	FS102-03
	6mL	FS102-03B
	3mL 50mL Pectinase	FS102-03K
	6mL & 50mL Pectinase – 1unit	FS102-03KB-200mg
	3mL & HPLC column & standard solution	FS102-03KLCS
	6mL & HPLC column & standard solution	FS102-03BKLCSP-200mg
	Pectinase enzymatic Solution 50mL – 1 unit	REA-001-50mL
	Patulin standard solution 1mL vial 100µg/mL	REA-PAT-1mL

**HPLC Column** : SilactHPLC LC-Patulin - 150x 2.1mm for patulin analysis ref: LC-Pat-150.2.1



## Analytes

- Ochratoxin A

## Matrices

- Wheat, Maize, red and white Wine, Several spices (Paprika, Pepper, ginger...), Coco, Human urine...

## Advantages

- Storage at RT like all AFFINIMIP®SPE
- Simple protocols
- Fluorescence or LC-MS detection

## Regulation

**Europe** - 1881/2006/EC: foodstuffs : 5µg/kg in raw cereal grains, 30µg/kg in spices, 2µg/kg in wine....  
**Codex Alimentarius** - (193-1995): 5µg/Kg for wheat, barley, rye...

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
AFFINIMIP® SPE Ochratoxin A	3mL	FS101-03
	6mL	FS101-03B
	10mL LRC	FS101-03LRC
	3mL for automate GERSTEL	FS101-03-GER
	3mL for automate GILSON	FS101-03-GIL
	6mL for automate GERSTEL	FS101-03B-GER
	6mL for automate GILSON	FS101-03B-GIL

*Visit our website to know more about this product and its formats.  
 Custom designed formats are available on demand.*



## Analytes

- Deoxynivalenol (*a.k.a* Vomitoxin)
- 3-AcetylDON
- 15-AcetylDON

## Matrices

- Oat, wheat, corn, baby food, meat, animal feed...

## Advantages

- UV and LC-MS detection

## Regulation

**Europe** - 1126/2007/EC : 1750µg/kg for unprocessed maize, 200µg/kg for babyfood cereals.

**US FDA** : 1 ppm in finished wheat products for human consumption

**Codex Alimentarius** - (193-1995): 1000µg/kg for flour maize, wheat, barley, 200µg/kg for babyfood cereals.

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
AFFINIMIP® SPE Deoxynivalenol	6mL -100mg in food and babyfood	FS117-03B
	6mL – 200mg in feed	FS117-03B-200mg

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*



## Analytes

- Zearalenone, Zearalanone,  $\alpha$ -zearalenol,  $\alpha$ -zearalanol,  $\beta$ -zearalenol,  $\beta$ -zearalanol

## Matrices

- Wheat, Maize, Cereal-based baby food, Edible corn oil and Rice

## Advantages

- Fluorescence or LC-MS detection

## Regulation

**Europe** - 1881/2006/EC : foodstuffs : 100 $\mu$ g/kg in cereals and 20 $\mu$ g/kg in maize-based babyfood.

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE</b> <b>Zearalenone</b>	3mL	FS100-03
	10mL LRC	FS100-03LRC
	3mL for automate GERSTEL	FS100-03-GER
	3mL for ASPEC automate GILSON	FS100-03-GIL

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*



## Analytes

- Glyphosate
- AMPA
- Glufosinate

## Matrices

- Waters: Geothermal, mineral, river or underground
- Foodstuffs: cereals, honey, tea, juices, cannabis

## Advantages

- NO DERIVATIZATION** required to extract the analytes
  - Fluorescence (with derivatization), LC-MS or Capillary Electrophoresis - UV detection
- Fast and simple protocol

## Regulation

In Europe, Glyphosate is a debatable active substance, it has been thoroughly assessed by Member States and the European Food Safety Authority (EFSA). In December 2017: the Commission has adopted an act to renew the approval of glyphosate for 5 years

**Codex alimentarius** : 50µg/kg in meat or milk

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE</b> <b>Glyphosate</b>	3mL	FS113-03.IP
	6mL	FS113-03B
	12mL	FS113-03C

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*

# AFFINIMIP® SPE Picolinic Herbicides



## Analytes

- Picloram
- Aminopyralid
- Clopyralid

## Matrices

- Water, Compost, Cereal, Soil...

## Advantages

- Fast, short and easy protocol

## Regulation

**Europe** - 322/2012/EC: 0.1 to 5 mg Clopyralid/kg in cereals  
 2017/171/EC: 10µg Aminopyralid/kg in nuts, vegetable, oilseeds...

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP®SPE Picolinic Herbicides</b>	3mL	FS115-03
	6mL	FS115-03B
	10mL LRC	FS115-03LRC
	12mL	FS115-03C

*Visit our website to know more about this product and its formats.  
 Custom designed formats are available on demand.*



## Analytes

- Benzo[a]anthracen B[a]A;  
Benzo[a]pyren B[a]P; Benzo[a]  
fluoranthen B[a]F; Chrysen (CHR),  
etc.

## Matrices

- Edible oils (colza, oliva...), fatty food

## Advantages

- LC-MS, HPLC/UV, Fluorescence
- Fast and simple protocol

## Regulation

**Europe** (EC 835/2011) 2 ng/g benzo[a]pyrene individually, and 10 ng/g benzo[a]pyrene, benzo[b]fluoranthene, chrysene and benzo[a]anthracene combined for edible oils

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE PAHs</b>	3mL	FS119-03-NG
	6mL	FS119-03B-NG
	10mL LRC	FS119-03LRC-NG
	12mL	FS119-03C-NG
	96 well plate – 1 unit	FS119-1.96W

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*

# AFFINIMIP® SPE Tetracyclines



## Analytes

- Tetracycline
- Chlortetracycline
- Oxytetracycline, their epimers
- Doxycycline.

## Matrices

- Meat, Tissues, Animal source foods, milk

## Advantages

- UV detection

## Regulation

**Europe** : 37/2010/EU  
100µg/kg in muscle or milk  
**Codex alimentarius** : 100  
µg/kg of milk; 200µg/kg of  
muscle

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE Tetracyclines</b>	1mL	FS112-03A
	3mL	FS112-03
	LRC 10mL	FS112-03LRC
	96 well plate – 1 unit	FS112-1.96W (1/pk)
	Multipurpose sampler automate (GERSTEL)-1mL	FS112-03A-GER
	ASPEC automate (GILSON)-1mL	FS112-03A-GIL
	ASPEC automate (GILSON)-3mL	FS112-03-GIL

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*

# AFFINIMIP® SPE Chloramphenicol



## Analytes

- Chloramphenicol

## Matrices

- Honey, Milk, Shrimp, Bovine Urine

## Advantages

- Very low limit of detection
- LC-MS

## Regulation

**Europe** - 2003/181/EC:  
minimum required  
performance 0.3µg/kg in  
residue of animal origin  
**US FDA** : prohibited

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE</b> <b>Chloramphenicol</b>	1mL	FS110-03A
	3mL	FS110-03
	LRC 10mL	FS110-03LRC
	96 well plate – 1 unit	FS110-1.96W (1/pk)

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*

# AFFINIMIP® SPE Aminoglycosides



## Analytes

- Spectinomycin, Hygromycin B, Streptomycin, Dihydrostreptomycin, Amikacin, Kanamycin A, Apramycin, Paromomycin, Tobramycin, Sisomicin, Gentamicin C1a, Gentamicin C2, Neomycin B

## Matrices

- Tissus muscles, milk, fish, egg

## Advantages

- **SIMULTANEOUS** analysis by LC-MS detection

## Regulation

<b>Europe</b> – 37/2010/EC:	50µg
total Gentamycins/kg	100µg
Kanamycin A/kg,	300µg
Spectinomycin/kg,	500µg
Dihydrostreptomycin in muscle	
<b>Codex Alimentarius</b> -	100µg
total Gentamycins/kg,	500µg
Spectinomycin/kg,	600µg
Dihydrostreptomycin	+
Streptomycin in muscle	

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE</b> <b>Aminoglycosides</b>	3 mL	FS124-03
	6mL	FS124-03B

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*

# AFFINIMIP® SPE Estrogens



## Analytes

- Broad family of natural and synthetic estrogens

## Matrices

- Waters, river water and sediment, Plasma, treated sewage, animal body fluid

## Advantages

- LC-MS, GC/MS

## Regulation

**Europe** : COM(2011)876 proposed Ethinylestradiol and 17  $\alpha$ -Ethinylestradiol as priority substances in water

**USA:** 17  $\alpha$ -Ethinylestradiol, Estriol, Estrone and Ethinylestradiol in Contaminant candidate list CCL3

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE Estrogens</b>	1mL	FS104-03A
	3mL	FS104-03
	6mL	FS104-03B
	10mL LRC	FS104-03LRC
	12mL	FS104-03C
	96 well plate – 1 unit	FS104-1.96W (1/pk)
	Reversible cart. – 0.7mL	FS104-03Rev1
	Reversible cart. – 2 mL	FS104-03Rev2

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*



## Analytes

- Bisphenols such as Bisphenol A and closely 18 related structures

## Matrices

- Water, milk (infant formula), powdered infant formula, canned food, vegetable puree for infant, Beer, urine, ...

## Advantages

- SIMULTANEOUS analysis of bisphenols analogs
- Broad range of solid and liquid foods tested
- LC-MS, GC/MS, fluorescence

## Regulation

**Europe 2011/8/EU** : forbidden in infant feeding bottle  
 Specific migration limit of 0.6mg/kg in food from packaging  
 Forbidden in food materials in France

Designation	Description	50 c/box
<b>AFFINIMIP® SPE Bisphenols</b>	3mL (PP)	FS106-03
	6mL (PP)	FS106-03B
	6mL (Glass)	FS106-03G
	10mL LRC (PP)	FS106-03LRC
	12mL (PP)	FS106-03C
	96 well plate – 1 unit	FS106-1.96W (1/pk)
	Multipurpose sampler automate (GERSTEL) -3mL	FS104-03-GER
	ASPEC automate (GILSON) - 3mL	FS104-03-GIL

HPLC Column : Silact LC-Bisphenol A - 150x 2.1mm for BPA analysis  
 ref: LC-BPA-150.2.1

*Visit our website to know more about this product and its formats.  
 Custom designed formats are available on demand.*



## Analytes

- Zeranol
- Zearalenone
- β-Zearalanol (Taleranol)
- α and β-Zearalenol
- Zearalanone
- Resorcylic acid lactones

## Matrices

- Meat, Urine, Tissues, Plasma

## Advantages

- LC-MS, HPLC-Fluorescence

## Regulation

Its use is banned in several countries (e.g. European directive 96/22/EC and in China).

**Codex alimentarius:** 2µg/kg in muscle

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
AFFINIMIP® SPE Zeranol Residues	3mL	FS105-03
	LRC 10mL	FS105-03LRC

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*



## Analytes

- Parabens, carnosic acid, hydroxylated PAHs
- Tocopherols
- Nitrophenols
- Chlorophenols
- Catechins...

## Matrices

- Food, cosmetic, wine, meat

## Advantages

- Extraction of broad range of phenolics compounds
- LC-MS, HPLC/UV

Designation	Description	50 c/box
<b>AFFINIMIP® SPE Phenolics</b>	3mL	FS103-03
	6mL	FS103-03B
	10mL LRC	FS103-03LRC
	12mL	FS103-03C
	96 well plate – 1 unit	FS103-1.96W (1/pk)
	Reversible cart. – 0.7mL	FS103-03Rev1
	Reversible cart. – 2 mL	FS103-03Rev2

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*



## Analytes

- **Total and free NNAL** (4-(methyl nitrosamino) -1-(3-pyridyl) -1-butanol)

## Matrices

- Urine

## Advantages

- LC-MS & LC-MS/MS

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE NNAL</b>	3mL	DG103-03
	10mL LRC	DG103-03LRC
	96 well plate – 1 unit	DG103-1.96W (1/pk)

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*



## Analytes

- Amphetamine
- Methamphetamine derivatives

## Matrices

- Urine, Serum

## Advantages

- LC-MS

## Regulation

Several US and European states have set up a cut-off value in urine or blood [e.g. France and Virginia (respectively 50ng/mL and 100ng/mL of blood)].

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE</b> <b>Amphetamines</b>	3mL	DG102-03
	10mL LRC	DG102-03LRC
	96 well plate – 1 unit	DG102-1.96W (1/pk)

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*



## Analytes

- Dopamine (DA)
- Norepinephrine or Noradrenaline (NA)
- Epinephrine or Adrenaline (A)

## Matrices

- Plasma , Serum

## Advantages

- LC-MS, HPLC/UV

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE</b> <b>Catecholamines</b>	1mL	DG100-03A
	3mL	DG100-03
	10mL LRC	DG100-03LRC
	12mL	DG100-03C
	96 well plate – 1 unit	DG100-1.96W (1/pk)

*Visit our website to know more about this product and its formats.  
 Custom designed formats are available on demand.*

# AFFINIMIP® SPE Metanephrines



## Analytes

- Metanephrine
- Normetanephrine
- 3-methoxytyramine

## Matrices

- Plasma , Serum

## Advantages

- LC-MS, HPLC/UV

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
<b>AFFINIMIP® SPE Metanephrines</b>	1mL	DG101-03A
	3mL	DG101-03
	10mL LRC	DG101-03LRC
	96 well plate – 1 unit	DG101-1.96W (1/pk)

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*



## Analytes

- Tamoxifen
- its metabolite: 4-HydroxyTamoxifen

## Matrices

- Biofluids such as urine

## Advantages

- LC-MS, HPLC/UV

## Regulation

Since January 2000, Tamoxifen has been included in the list of prohibited substances by the International Olympic Committee.

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
AFFINIMIP® SPE Tamoxifen	3mL	PH101-03

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*



## Analytes

- Tetrahydrocannabinol (THC) and its metabolites
- THC-OH
- THC-COOH

## Matrices

- Biofluids such as urine, saliva, blood

## Advantages

- LC-MS

## Regulation

This drug is forbidden in most countries in the world. However, some countries or US states make possible the use of this product for medical use (e.g. Arizona, Florida, Louisiana...) and for even less countries/states for recreational use (e.g. Canada, California, Alaska, Massachusetts...).

Designation	Description	50 c/box
AFFINIMIP® SPE Cannabis	3mL	PH123-03
	6mL	PH123-03B
	10mL LRC	PH123-03LRC

*Visit our website to know more about this product and its formats.  
Custom designed formats are available on demand.*

# AttractSPE® PFAS FOR PERFLUORINATED COMPOUNDS

AttractSPE® PFAS specially developed and dedicated to determine perfluorinated compounds.



## Analytes

- Perfluorobutanoic acid (PFBA), Perfluoropentanoic acid (PFPeA), Perfluorohexanoic acid (PFHxA), Perfluoroheptanoic acid (PFHpA), Perfluorooctanoic acid (PFOA), Perfluorononanoic acid (PFNA), Perfluorodecanoic acid (PFDA), Perfluorotetradecanoic acid (PFTA), Perfluorobutanesulfonic acid (PFBS), Perfluorohexane sulfonic acid (PFHxS), Perfluorooctanesulfonic acid (PFOS)

## Matrices

Water

## Advantages

- LC-MS/MC

*Application notes and publications are available on our website.*

Designation	Description	50 c/box
AttractSPE® PFAS	6mL – 150mg	PFAS-50.S.6.150
	6mL – 200mg	PFAS-50.S.6.200
	6mL – 500mg -80µm	PFAS-50.S.6.500GP

### HPLC columns for perfluorinated compounds (PFAS):

SilactHPLC DELAY - PFAS column: Delay column for analysis of PFAS - 50x2.1mm - ref : **DELAY-PFAS-50.2.1**

SilactHPLC LC - PFAS column- 150x2.1mm - ref : **LC-PFAS-150.2.1**

*Visit our website to know more about this product and its formats.*

*Custom designed formats are available on demand.*

# AttractSPE® NON-TARGET SCREENING & AttractSPE® CHELATING FOR METAL RETENTION

**AttractSPE® NON-TARGET SCREENING** : a mixture of sorbents for Non-target screening

**AttractSPE® NON-TARGET SCREENING** are combination of sorbents to detect the presence of not regulated harmful substances (including degradation by-products and metabolites)

**AttractSPE® CHELATING** : A chelating resin for the extraction of metal ions

**AttractSPE® CHELATING** is a polymer resin containing functional groups which particularly suits for the extraction of transition metal ions (Hg<sup>2+</sup>, Cu<sup>2+</sup>, Pb<sup>2+</sup>, Fe<sup>2+</sup>...) and Alkaline earth metals ions in water even in highly concentrated salt solution. A high selectivity towards metal ions is obtained by varying the pH.

Cartridges format, Description	AttractSPE® Non-Target Screening (50/pk)
6mL , 300mg Mixture of HLB-WCX-WAX	Screening1-50.S.6.300
6mL , 300mg Mixture of HLB-SCX-SAX	Screening2-50.S.6.300
6mL , 300mg Mixture of PSDVB-WCX-WAX	Screening3-50.S.6.300
6mL , 300mg Mixture of PSDVB-SCX-SAX	Screening4-50.S.6.300
47mm or 90mm Disks	See page <b>AttractSPE® Disks for environmental applications</b>

Cartridges format, Sorbent amount	AttractSPE® CHELATING (50/pk)
1mL , 30mg	MET-100.S.1.30 (100/pk)
3mL, 60mg	MET50.S.3.60
6mL, 100mg	MET-50.S.6.200
6mL, 500mg	MET-50.S.6.500
96 wells Plate	MET-1.96W.30 (1/pk)
Reversible 0.7mL, 30mg	MET-50.REV.1.N10
Reversible 0.7mL, 100mg	MET-50.REV.1.F
Reversible 2mL, 800mg	MET-50.REV.2.F
47mm or 90mm Disks	See page <b>AttractSPE® Disks for environmental applications</b>

**SilactSPE DAU** particularly suits for basic drugs of abuse determination (Amphetamines, opioids...) in complex matrices such as human urine.

## Product Information

**Silica type** : 60 Å, 500 m<sup>2</sup>/g, 40-63 µm

Cartridges format, Sorbent amount	#/box	SilactSPE DAU
1mL, 50mg	100	DAU-100.S.1.50
1mL, 100mg	100	DAU-100.S.1.100
3mL, 200mg	50	DAU-50.S.3.200
3mL, 500mg	50	DAU-50.S.3.500
6mL, 500mg	50	DAU-50.S.6.500
6mL, 1g	50	DAU-50.S.6.1g
10mL LRC, 500mg	50	DAU-50.LRC.10.500
12mL, 2g	20	DAU-20.S.12.2g

**SilactSPE Organotins** : for organotins analysis following ISO 17353 and ISO 23161 methods

**SilactSPE Organotins** is a column for the analysis of these compounds in water, soil, sediment, sludge and waste (soil-like material).

This column can be used for the analysis of:

- **Organotins cations** (Monobutyltin BuSn<sup>3+</sup>; Dibutyltin Bu<sub>4</sub>Sn<sup>3+</sup>; Tributyltin Bu<sub>3</sub>Sn<sup>+</sup>; Monooctyltin OcSn<sup>3+</sup>; Dioctyltin Oc<sub>2</sub>Sn<sup>2+</sup>; Triphenyltin Ph<sub>3</sub>Sn<sup>+</sup>; Tricyclohexyltin Cy<sub>3</sub>Sn<sup>+</sup>)

- **Peralkylated organotin** (Tetrabutyltin Bu<sub>4</sub>Sn)

- **Methyltin compounds** (Monomethyltin MeSn<sup>3+</sup>; Dimethyltin Me<sub>2</sub>Sn<sup>2+</sup>; Trimethyltin Me<sub>3</sub>Sn<sup>+</sup>)

Cartridges format, Sorbent amount	#/box	SilactSPE Organotins
25mL, 3g Na <sub>2</sub> SO <sub>4</sub> + 5g Silica (3% water)	25	OSn-25.S.25.3g.5g
	50	Osn-50.S.25.3g.5g

# SPE FOR POLYCYCLIC AROMATIC HYDROCARBONS (PAHs)

## AFFINIMIP®SPE

### PAHs

For the cleanup of PAHs in **FATTY food and liquid** such as oil

Molecularly imprinted polymer for PAHs.

## SilactSPE

### CN/SiOH

For the cleanup of PAHs in **SOIL**

A two layer sorbents with cyano modified silica and silica sorbents

## AttractSPE®

### HLB

For the cleanup of PAHs in **WATER**

HLB

Product	Vol	Sorbent	25 cartridges/box	50 cartridges/box
<b>SilactSPE</b> <b>CN/SiOH</b>	3mL	500mg/1g	CNSiOH-25.S.3.500.1g	CNSiOH-50.S.3.500.1g
	6mL	500mg/1g	CNSiOH-25.S.6.500.1g	CNSiOH-50.S.6.500.1g
	6mL glass	500mg/1g	CNSiOH-25.G.6.500.1g	CNSiOH-50.G.6.500.1g
<b>AFFINIMIP®</b> <b>SPE PAHs</b>	3mL			FS119-03-NG
<b>AttractSPE®</b> <b>HLB</b>	6mL	200mg		HLB-50.S.6.200

# SPE FOR INTERFERENCES REMOVAL

## AttractSPE®SAX-HCO3

For the removal of anionic contaminants and neutralization of acidic samples

**AttractSPE®SAX-HCO3** is a strong anion exchange sorbent with hydrogenocarbonate anion as counterion. It is used for the removal of anionic contaminants from sample matrices and for the neutralization of highly acidic samples.

### Product Information

**PS-DVB type:** 40 µm, 60 Å, 600 m<sup>2</sup>/g, 0.3 meq/g

## AttractSPE®PS-H

For the removal of alkaline earth and transition metals ions and to neutralize basic samples.

**AttractSPE®PS-H** is a strong cation exchange sorbent in the H form. It is used for the removal of alkaline earth and transition metals ions and to neutralize basic samples.

### Product Information

**PS-DVB polymer type:** 60 Å, 600m<sup>2</sup>/g, 1meq/g, 40 µm

Cartridges format, Sorbent amount	# /box	AttractSPE®SAX-HCO3	AttractSPE®PS-H
1mL	100	SAX-HCO3-100.S.1.30	
3mL, 60mg	50	SAX-HCO3-50.S.3.60	PSH-50.S.3.60
6mL, 200mg	50	SAX-HCO3-50.S.6.200	PSH-50.S.6.200
6mL, 500mg	50	SAX-HCO3-50.S.6.500	PSH-50.S.6.500
96 well Plate	1	SAX-HCO3-1.96W.30	PSH-1.96W.30
Reversible 0.7mL, 100mg	25	SAX-HCO3-25.REV.1.F	PSH-25.REV.1.F
	50	SAX-HCO3-50.REV.1.F	PSH-50.REV.1.F
Reversible 2mL, 800mg	25	SAX-HCO3-25.REV.2.F	PSH-25.REV.2.F
	50	SAX-HCO3-50.REV.2.F	PSH-50.REV.2.F

## SPE FOR INTERFERENCES REMOVAL

### AttractSPE® PS-Ag

Removal of halide ions (**chloride, bromide, and iodide**) by precipitation

**Strong cation exchange sorbent with silver cation** as counterion.

### AttractSPE® PS-Ba

Removal of sulfate ions by precipitation

**Strong cation exchange sorbent with barium cation** as counterion

### SilactSPE HydroxyApatite

Removal of **chloride, fluoride, lanthanide & carbonate ions**


Hydroxyapatite is a mineral compound of structure  $\text{Ca}_{10}(\text{PO}_4)_6(\text{OH})_2$

Cartridges format, Sorbent amount	# /box	AttractSPE® PS-Ag	AttractSPE® PS-Ba	SilactSPE HydroxyApatite
1mL, 30mg (50mg for HAp)	100	PSAg-100.S.1.30	PSBa-100.S.1.30	HAp-100.S.1.50
3mL, 60mg (200mg for HAp)	50	PSAg-50.S.3.60	PSBa-50.S.3.60	HAp-50.S.3.200
6mL, 200mg	50	PSAg-50.S.6.200	PSBa-50.S.6.200	-
6mL, 500mg	50	PSAg-50.S.6.500	PSBa-50.S.6.500	HAp-50.S.6.500
Reversible 0.7mL, 400mg	25	PSAg-25.REV.1.F	PSBa-25.REV.1.F	
	50	PSAg-50.REV.1.F	PSBa-50.REV.1.F	HAp-50.REV.1.F



# AttractSPE® Cartridges

**A broad range of chemistry and format for your application**




## POLYMERIC - BASED AttractSPE® CARTRIDGES

**AttractSPE®** are based on polymeric sorbents dedicated to the extraction of compounds from complex matrices. **AttractSPE®** cartridges provide the advantages of robustness, simplified method development, wide applicability and are not affected by drying out. The combination of the water-wettable optimised surface chemistry, high surface area and pH stability ensures high reproducible recoveries for a wide range of analytes. Two particle sizes (40 and 60µm) are available for a more flexible implementation of these products.

The choice of the suitable **AttractSPE®** columns can be done by using the following method:

- Check if a method already exists on our application notebook
- Determination of the nature of the analyte (neutral, acid, base)
- Determination of the pKa
- Choice of the **AttractSPE®** columns by using the following chart and application of the general protocol described on the instruction sheet
- Analysis of the recovery yields





**AttractSPE<sup>®</sup> HLB** is an **uncharged Hydrophilic and Lipophilic sorbent** interacting with both, hydrophilic and hydrophobic interactions. It particularly suits for the extraction of a wide range of analytes (polar, apolar, neutral, acid, basic...)

### Product Information

**Particle diameter range** : 40 and 80µm

**Pore size**: 70 Å

**Surface area**: 800 m<sup>2</sup>/g

**Storage** : Ambient temperature

Format, amount	#/box	AttractSPE <sup>®</sup> HLB -40µm	AttractSPE <sup>®</sup> HLB -80µm
1mL, 10mg, 100/pk	100	HLB-100.S.1.10	
1mL, 30mg	100	HLB-100.S.1.30	HLB-100.S.1.30GP
3mL, 60mg	50	HLB-50.S.3.60	HLB-50.S.3.60GP
3mL, 100mg	50	HLB-50.S.3.100	HLB-50.S.3.100GP
6mL, 150mg	50	HLB-50.S.6.150	HLB-50.S.6.150GP
6mL, 200mg	50	HLB-50.S.6.200	HLB-50.S.6.200GP
6mL, 500mg	50	HLB-50.S.6.500	HLB-50.S.6.500GP
10mL LRC, 60mg	50	HLB-50.LRC.10.60	
12mL, 500mg	25	HLB-25.S.12.500	
12mL, 1g	25	HLB-25.S.12.1000	
12mL, 2g	25		HLB-25.S.12.2000GP
20mL, 1g	25	HLB-25.S.20.1000	
20mL, 2g	25		HLB-25.S.20.2000
96 well Plate, 10mg	1	HLB-1.96W.10	
96 well Plate, 30mg	1	HLB-1.96W.30	HLB-1.96W.30GP
96 well Plate, 60mg	1	HLB-1.96W.60	HLB-1.96W.60GP
Reversible 0.7mL, 30mg	50	HLB-50.REV.1.N10	
Reversible 0.7mL, 100mg	50	HLB-50.REV.1.F	
Reversible 2mL, 225mg	50	HLB-50.REV.2.N10	
Disks cartridges (sorbent is a SPE disk)		See page <b>AttractSPE<sup>®</sup> Disks Cartridges</b> High capacity for a minimal elution volume	

## Mixed-mode SPE for extraction of strong acid analytes

**AttractSPE® WAX** is a weak anion exchange sorbent interacting with the analytes via a mixed mode retention mechanism, ion exchange with weak basic functional groups and reversed phase. It particularly suits for the extraction of strong acids.

## Product Information

**Particle diameter range:** 40 and 80µm

**Pore size:** 60 Å

**Surface area:** 650 m<sup>2</sup>/g

**Ionic capacity :** 0.5 meq/g

**Storage:** Ambient temperature

Format, amount	#/box	AttractSPE®WAX-40µm	AttractSPE® WAX -80µm
1mL, 10mg, 100/pk	100	WAX-100.S.1.10	
1mL, 30mg	100	WAX-100.S.1.30	WAX-100.S.1.30GP
3mL, 60mg	50	WAX-50.S.3.60	WAX-50.S.3.60GP
3mL, 100mg	50	WAX-50.S.3.100	WAX-50.S.3.100GP
6mL, 150mg	50	WAX-50.S.6.150	WAX-50.S.6.150GP
6mL, 200mg	50	WAX-50.S.6.200	WAX-50.S.6.200GP
6mL, 500mg	50	WAX-50.S.6.500	WAX-50.S.6.500GP
10mL LRC, 60mg	50	WAX-50.LRC.10.60	
12mL, 500mg	25	WAX-25.S.12.500	
12mL, 1g	25	WAX-25.S.12.1000	
12mL, 2g	25		WAX-25.S.12.2000GP
20mL, 1g	25	WAX-25.S.20.1000	
20mL, 2g	25		WAX-25.S.20.2000
96 well Plate, 10mg	1	WAX-1.96W.10	
96 well Plate, 30mg	1	WAX-1.96W.30	WAX-1.96W.30GP
96 well Plate, 60mg	1	WAX-1.96W.60	WAX-1.96W.60GP
Reversible 0.7mL, 30mg	50	WAX-50.REV.1.N10	
Reversible 0.7mL, 100mg	50	WAX-50.REV.1.F	
Reversible 2mL, 225mg	50	WAX-50.REV.2.N10	

**For perfluorinated compounds, please use **AttractSPE® PFAS****

## Mixed-mode SPE for extraction of strong basic analytes

AttractSPE® WCX is a weak cation exchange sorbent interacting with the analytes via a mixed mode retention mechanism, ion exchange with weak acid functional groups and reversed phase. It particularly suits for the extraction of strong bases and quaternary amines.

## Product Information

**Particle diameter range :** 40 and 80µm

**Pore size:** 70 Å

**Surface area:** 850 m<sup>2</sup>/g

**Ionic capacity:** 0.77meq/g

**Storage :** Ambient temperature

Format, amount	#/box	AttractSPE®WCX-40µm	AttractSPE® WCX -80µm
1mL, 10mg, 100/pk	100	WCX-100.S.1.10	
1mL, 30mg	100	WCX-100.S.1.30	WCX-100.S.1.30GP
3mL, 60mg	50	WCX-50.S.3.60	WCX-50.S.3.60GP
3mL, 100mg	50	WCX-50.S.3.100	WCX-50.S.3.100GP
6mL, 150mg	50	WCX-50.S.6.150	WCX-50.S.6.150GP
6mL, 200mg	50	WCX-50.S.6.200	WCX-50.S.6.200GP
6mL, 500mg	50	WCX-50.S.6.500	WCX-50.S.6.500GP
10mL LRC, 60mg	50	WCX-50.LRC.10.60	
12mL, 500mg	25	WCX-25.S.12.500	
12mL, 1g	25	WCX-25.S.12.1000	
12mL, 2g	25		WCX-25.S.12.2000GP
20mL, 1g	25	WCX-25.S.20.1000	
20mL, 2g	25		WCX-25.S.20.2000
96 well Plate, 10mg	1	WCX-1.96W.10	
96 well Plate, 30mg	1	WCX-1.96W.30	WCX-1.96W.30GP
96 well Plate, 60mg	1	WCX-1.96W.60	WCX-1.96W.60GP
Reversible 0.7mL, 30mg	50	WCX-50.REV.1.N10	
Reversible 0.7mL, 100mg	50	WCX-50.REV.1.F	
Reversible 2mL, 225mg	50	WCX-50.REV.2.N10	

## Mixed-mode SPE for extraction of weak acid analytes

**AttractSPE® SAX** is a strong anion exchange sorbent interacting with the analytes via a mixed mode retention mechanism, ion exchange with strong basic functional groups and reversed phase. It particularly suits for the extraction of weak acids.

## Product Information

**Particle diameter range** : 40 and 80µm

**Pore size**: 60 Å

**Surface area**: 600 m<sup>2</sup>/g

**Ionic capacity**: 0.3 meq/g

**Storage** : Ambient temperature

Format, amount	#/box	AttractSPE® SAX -40µm	AttractSPE® SAX -80µm
1mL, 10mg, 100/pk	100	SAX-100.S.1.10	
1mL, 30mg	100	SAX-100.S.1.30	SAX-100.S.1.30GP
3mL, 60mg	50	SAX-50.S.3.60	SAX-50.S.3.60GP
3mL, 100mg	50	SAX-50.S.3.100	SAX-50.S.3.100GP
6mL, 150mg	50	SAX-50.S.6.150	SAX-50.S.6.150GP
6mL, 200mg	50	SAX-50.S.6.200	SAX-50.S.6.200GP
6mL, 500mg	50	SAX-50.S.6.500	SAX-50.S.6.500GP
10mL LRC, 60mg	50	SAX-50.LRC.10.60	
12mL, 500mg	25	SAX-25.S.12.500	
12mL, 1g	25	SAX-25.S.12.1000	
12mL, 2g	25		SAX-25.S.12.2000GP
20mL, 1g	25	SAX-25.S.20.1000	
20mL, 2g	25		SAX-25.S.20.2000
96 well Plate, 10mg	1	SAX-1.96W.10	
96 well Plate, 30mg	1	SAX-1.96W.30	SAX-1.96W.30GP
96 well Plate, 60mg	1	SAX-1.96W.60	SAX-1.96W.60GP
Reversible 0.7mL, 30mg	50	SAX-50.REV.1.N10	
Reversible 0.7mL, 100mg	50	SAX-50.REV.1.F	
Reversible 2mL, 225mg	50	SAX-50.REV.2.N10	
Disks cartridges (sorbent is a SPE disk)		See page <b>AttractSPE® Disks Cartridges</b> High capacity for a minimal elution volume	

## Mixed-mode SPE for extraction of weak basic analytes

**AttractSPE® SCX** is a strong cation exchange sorbent interacting with the analytes via a mixed mode retention mechanism, ion exchange with strong acid functional groups and reversed phase. It particularly suits for the extraction of weak bases.

## Product Information

**Particle diameter range** : 40 and 80 µm

**Pore size**: 60 Å

**Surface area**: 600 m<sup>2</sup>/g

**Ionic capacity**: 1meq/g

**Storage** : Ambient temperature

Format, amount	#/box	AttractSPE® SCX -40µm	AttractSPE® SCX -80µm
1mL, 10mg, 100/pk	100	SCX-100.S.1.10	
1mL, 30mg	100	SCX-100.S.1.30	SCX-100.S.1.30GP
3mL, 60mg	50	SCX-50.S.3.60	SCX-50.S.3.60GP
3mL, 100mg	50	SCX-50.S.3.100	SCX-50.S.3.100GP
6mL, 150mg	50	SCX-50.S.6.150	SCX-50.S.6.150GP
6mL, 200mg	50	SCX-50.S.6.200	SCX-50.S.6.200GP
6mL, 500mg	50	SCX-50.S.6.500	SCX-50.S.6.500GP
10mL LRC, 60mg	50	SCX-50.LRC.10.60	
12mL, 500mg	25	SCX-25.S.12.500	
12mL, 1g	25	SCX-25.S.12.1000	
12mL, 2g	25		SCX-25.S.12.2000GP
20mL, 1g	25	SCX-25.S.20.1000	
20mL, 2g	25		SCX-25.S.20.2000
96 well Plate, 10mg	1	SCX-1.96W.10	
96 well Plate, 30mg	1	SCX-1.96W.30	SCX-1.96W.30GP
96 well Plate, 60mg	1	SCX-1.96W.60	SCX-1.96W.60GP
Reversible 0.7mL, 30mg	50	SCX-50.REV.1.N10	
Reversible 0.7mL, 100mg	50	SCX-50.REV.1.F	
Reversible 2mL, 225mg	50	SCX-50.REV.2.N10	
Disks cartridges (sorbent is a SPE disk)		See page <b>AttractSPE® Disks Cartridges</b> High capacity for a minimal elution volume	

## Reversed phase SPE for extraction of hydrophobic analytes

**AttractSPE® DVB** is a polystyrene-divinylbenzene copolymer presenting a high hydrophobicity used as a reversed-phase. It particularly suits for the extraction of hydrophobic analytes.

## Product Information

**Particle diameter range** : 40

and 80µm

**Pore size**: 60 Å

**Surface area**: 600 m<sup>2</sup>/g

**Storage** : Ambient temperature

Format, amount	#/box	AttractSPE® DVB -40µm	AttractSPE® DVB -80µm
1mL, 10mg, 100/pk	100	DVB-100.S.1.10	
1mL, 30mg	100	DVB-100.S.1.30	DVB-100.S.1.30GP
3mL, 60mg	50	DVB-50.S.3.60	DVB-50.S.3.60GP
3mL, 100mg	50	DVB-50.S.3.100	DVB-50.S.3.100GP
6mL, 150mg	50	DVB-50.S.6.150	DVB-50.S.6.150GP
6mL, 200mg	50	DVB-50.S.6.200	DVB-50.S.6.200GP
6mL, 500mg	50	DVB-50.S.6.500	DVB-50.S.6.500GP
10mL LRC, 60mg	50	DVB-50.LRC.10.60	
12mL, 500mg	25	DVB-25.S.12.500	
12mL, 1g	25	DVB-25.S.12.1000	
12mL, 2g	25		DVB-25.S.12.2000GP
20mL, 1g	25	DVB-25.S.20.1000	
20mL, 2g	25		DVB-25.S.20.2000
96 well Plate, 10mg	1	DVB-1.96W.10	
96 well Plate, 30mg	1	DVB-1.96W.30	DVB-1.96W.30GP
96 well Plate, 60mg	1	DVB-1.96W.60	DVB-1.96W.60GP
Reversible 0.7mL, 30mg	50	DVB-50.REV.1.N10	
Reversible 0.7mL, 100mg	50	DVB-50.REV.1.F	
Reversible 2mL, 225mg	50	DVB-50.REV.2.N10	
Disks cartridges (sorbent is a SPE disk)		See page <b>AttractSPE® Disks Cartridges</b> High capacity for a minimal elution volume	



# INORGANIC-BASED SilactSPE CARTRIDGES

## A very large range of SPE sorbents

**SilactSPE** products are inorganic based sorbents SPE cartridges mainly alumina or modified silica.

**SilactSPE Silica or Alumina - based SPE cartridges** are silica- or alumina based phases and offer a broad range of chemically modified silica or alumina. This chemistry goes from very polar sorbent (bare silica) to hydrophobic sorbent (end-capped saturated hydrocarbon modified silica) passing through intermediate polarity (for instance, amino modified silica).

**SilactSPE** products are Silica-based and alumina-based sorbents available in different formats including SPE cartridges and 48- & 96-well plates, with different sorbents, and in bed weights up to 10 grams.

Reversed phase based silica	More polar Silica based phase		Normal phase
<b>C8</b> moderately hydrophobic	<b>SiWCX</b> Weak cation exchanger	<b>SiSCX</b> Strong cation exchanger	<b>Silica</b> Very polar
<b>Phenyl</b> moderately hydrophobic	<b>SiSAX</b> Strong anion exchanger	<b>Amino (SiWAX)</b> Weak anion exchanger	<b>Alumina (A, B, N)</b> Highly active
<b>C18</b> Strongly hydrophobic	<b>Cyano</b> Cyano propyl Polar phase	<b>PSA</b> primary secondary amine	<b>Florisil</b> polar –highly active – weakly basic



## Strongly hydrophobic and non-polar sorbent

It was recently developed as an innovative C18 phase characterized by a homogeneous coverage of the silane on the surface.

**SilactSPE C18** particularly suits for the extraction of acidic, neutral and basic compounds from aqueous solutions, various organic compounds from water, and drugs and metabolites from physiological fluids.


## Product Information

**Loading** : 15-17 % C

**Endcapping** : Yes for **SilactSPE C18**,

No for **SilactSPE C18NEC**

**Silica type**: 60 Å, 500 m<sup>2</sup>/g, 40-63 µm



**SilactSPE C18 end-capped**



**SilactSPE C18NEC not end-capped**

Cartridges format, Sorbent amount	#/box	SilactSPE C18 (end capped)	SilactSPE C18 NEC (not end capped)
1mL, 50mg	100	C18-100.S.1.50	C18nec-100.S.1.50
1mL, 100mg	100	C18-100.S.1.100	C18nec-100.S.1.100
3mL, 200mg	50	C18-50.S.3.200	C18nec-50.S.3.200
3mL, 500mg	50	C18-50.S.3.500	C18nec-50.S.3.500
6mL, 500mg	50	C18-50.S.6.500	C18nec-50.S.6.500
6mL, 1g	50	C18-50.S.6.1000	C18nec-50.S.6.1000
10mL LRC, 500mg	50	C18-50.LRC.10.500	C18nec-50.LRC.10.500
12mL, 2g	20	C18-20.S.12.2000	C18nec-20.S.12.2000
Reversible 0.7mL, 260mg	25	C18-25.REV.1.260	C18nec-25.REV.1.360
Reversible 2mL, 1g	25	C18-25.REV.2.1000	C18nec-25.REV.2.1000
96 well plate, 50mg	1	C18-1.96W.50	C18nec-1.96W.50
96 well plate, 100mg	1	C18-1.96W.100	C18nec-1.96W.100

## SilactSPE C8 & SilactSPE PHENYL

**SilactSPE C8:** Moderately hydrophobic and non-polar sorbent  
**Sorbent C8** is more selective than **Sorbent C18** for big compounds such as PAH, vitamin D, and oils as well as greasy compounds. It particularly suits for the extraction of extremely non-polar compounds.

### Product Information

**Loading :** 12 % C

**Endcapping :** Yes

**Silica type :** 60 Å, 500 m<sup>2</sup>/g, 40-63 µm

**SilactSPE Phenyl:** Moderately hydrophobic and non-polar sorbent  
 it particularly suits for the extraction of non-polar compounds with different selectivities through  $\pi$ - $\pi$  interactions including aromatic compounds and other non-polar phases.

### Product Information

**Loading :** 9 % C

**Endcapping :** Yes

**Silica type :** 60 Å, 500 m<sup>2</sup>/g, 40-63 µm

Cartridges format, Sorbent amount	#/box	SilactSPE C8	SilactSPE Phenyl
1mL, 50mg	100	C8-100.S.1.50	Phe-100.S.1.50
1mL, 100mg	100	C8-100.S.1.100	Phe-100.S.1.100
3mL, 200mg	50	C8-50.S.3.200	Phe-50.S.3.200
3mL, 500mg	50	C8-50.S.3.500	Phe-50.S.3.500
6mL, 500mg	50	C8-50.S.6.500	Phe-50.S.6.500
6mL, 1g	50	C8-50.S.6.1000	Phe-50.S.6.1000
10mL LRC, 500mg	50	C8-50.LRC.10.500	Phe-50.LRC.10.500
12mL, 2g	20	C8-20.S.12.2000	Phe-20.S.12.2000
Reversible 0.7mL, 260mg	25	C8-25.REV.1.260	Phe-25.REV.1.260
Reversible 2mL, 1g	25	C8-25.REV.2.1000	Phe-25.REV.2.1000
96 well plate, 50mg	1	C8-1.96W.50	Phe--1.96W.50
96 well plate, 100mg	1	C8-1.96W.100	Phe-1.96W.100

### SilactSPE Silica : Most polar sorbent

It presents a slightly acidic character and is used to extract various compounds from non-polar solvents through hydrogen bonding.

#### Product Information

**Silica type** : 60 Å, 500 m<sup>2</sup>/g, 40-63 µm

### SilactSPE Cyano : Moderately polar sorbent

It is used as a normal phase (less polar compared to silica). It particularly suits for the extraction of acidic, basic and neutral compounds from aqueous solutions. It is also used as a reversed-phase (less hydrophobic than C8 and C18).

#### Product Information

**Loading** : 7 % C

**Endcapping** : Yes

**Silica type** : 60 Å, 500 m<sup>2</sup>/g, 40-63 µm

Cartridges format, Sorbent amount	#/box	SilactSPE Silica	SilactSPE Cyano
1mL, 50mg	100	Si-100.S.1.50	CN-100.S.1.50
1mL, 100mg	100	Si-100.S.1.100	CN-100.S.1.100
3mL, 200mg	50	Si-50.S.3.200	CN-50.S.3.200
3mL, 500mg	50	Si-50.S.3.500	CN-50.S.3.500
6mL, 500mg	50	Si-50.S.6.500	CN-50.S.6.500
6mL, 1g	50	Si-50.S.6.1000	CN-50.S.6.1000
10mL LRC, 500mg	50	Si-50.LRC.10.500	CN-50.LRC.10.500
12mL, 2g	20	Si-20.S.12.2000	CN-20.S.12.2000
Reversible 0.7mL	25	Si-25.REV.1.240 for 240mg	CN-25.REV.1. 260 for 260mg
Reversible 2mL	25	Si-25.REV.2.900 for 900mg	CN-25.REV.2.1000 for 1000mg
96 well plate, 50mg	1	Si-1.96W.50	CN-1.96W.50
96 well plate, 100mg	1	Si-1.96W.100	CN-1.96W.100

## SilactSPE Amine (SiWAX) & SilactSPE PSA

**SilactSPE Amine (SiWAX):** Weak anion exchanger silica-based sorbent

**SilactSPE Amino** avoids irreversible retention of acidic molecules ( $pK_a < 3$ ) and particularly suits for the separation of peptides, drugs and metabolites from physiological fluids, poly- and monosaccharides and structural isomers.

### Product Information

**Loading :** 1.6 mmol/g

**Endcapping :** Yes

**Silica type :** 60 Å, 500 m<sup>2</sup>/g, 40-63 µm

**SilactSPE PSA:** Weak anion exchanger silica-based sorbent

Less polar sorbent than **SilactSPE Amine** used for its replacement with analytes that are too strongly retained on an amine phase.

### Product Information

**Endcapping :** Yes

**Silica type :** 60 Å, 500 m<sup>2</sup>/g, 40-63 µm

Cartridges format, Sorbent amount	#/box	SilactSPE Amine or SiWAX	SilactSPE PSA
1mL, 50mg	100	NH2-100.S.1.50	PSA-100.S.1.50
1mL, 100mg	100	NH2-100.S.1.100	PSA-100.S.1.100
3mL, 200mg	50	NH2-50.S.3.200	PSA-50.S.3.200
3mL, 500mg	50	NH2-50.S.3.500	PSA-50.S.3.500
6mL, 500mg	50	NH2-50.S.6.500	PSA-50.S.6.500
6mL, 1g	50	NH2-50.S.6.1000	PSA-50.S.6.1000
10mL LRC, 500mg	50	NH2-50.LRC.10.500	PSA-50.LRC.10.500
12mL, 2g	20	NH2-20.S.12.2000	PSA-20.S.12.2000
Reversible 0.7mL, 260mg	25	NH2-25.REV.1.260	PSA-25.REV.1.260
Reversible 2mL, 1000mg	25	NH2-25.REV.2.1000	PSA-25.REV.2.1000
96 well plate, 50mg	1	NH2-1.96W.50	PSA--1.96W.50
96 well plate, 100mg	1	NH2-1.96W.100	PSA-1.96W.100

**SilactSPE SiWCX:** Weak cation exchanger silica-based sorbent with carboxylic acid.

**SilactSPE SiWCX** particularly suits to extract strong basic molecules (pKa>9).

**Product Information**

**Loading :** 1.6 mmol/g

**Endcapping :** Yes

**Silica type :** 60 Å, 500 m<sup>2</sup>/g, 40-63 µm

**SilactSPE SiSCX:** Strong cation exchanger silica-based sorbent positively charged with tosic acid moieties.

**SilactSPE SiSCX** particularly suits to extract basic molecules (pKa 7-10)

**Product Information**

**Loading :** 0.8 mmol/g

**Endcapping :** Yes

**Silica type :** 60 Å, 500 m<sup>2</sup>/g, 40-63 µm

Cartridges format, Sorbent amount	#/box	SilactSPE SiWCX	SilactSPE SiSCX
1mL, 50mg	100	SiWCX-100.S.1.50	SiSCX-100.S.1.50
1mL, 100mg	100	SiWCX-100.S.1.100	SiSCX-100.S.1.100
3mL, 200mg	50	SiWCX-50.S.3.200	SiSCX-50.S.3.200
3mL, 500mg	50	SiWCX-50.S.3.500	SiSCX-50.S.3.500
6mL, 500mg	50	SiWCX-50.S.6.500	SiSCX-50.S.6.500
6mL, 1g	50	SiWCX-50.S.6.1000	SiSCX-50.S.6.1000
10mL LRC, 500mg	50	SiWCX-50.LRC.10.500	SiSCX-50.LRC.10.500
12mL, 2g	20	SiWCX-20.S.12.2000	SiSCX-20.S.12.2000
Reversible 0.7mL, 260mg	25	SiWCX-25.REV.1.260	SiSCX-25.REV.1.260
Reversible 2mL, 1000mg	25	SiWCX-25.REV.2.1000	SiSCX-25.REV.2.1000
96 well plate, 50mg	1	SiWCX-1.96W.50	SiSCX-1.96W.50
96 well plate, 100mg	1	SiWCX-1.96W.100	SiSCX-1.96W.100

**SilactSPE SiSAX:** Strong anion exchanger silica-based sorbent using trimethyl ammonium moieties.

**SilactSPE SiSAX** particularly suits to extract acidic molecules (pKa 3-5)

**Product Information**

**Loading :** 1.1 mmol/g

**Endcapping :** No

**Silica type :** 60 Å, 500 m<sup>2</sup>/g, 40-63 µm

**SilactSPE Carbonate**

**General base quencher**

**SilactSPE Carbonate** is the silica-bound equivalent of tetramethyl ammonium carbonate and is used as a general base to quench a reaction, free base amines in their ammonium salt form and to scavenge acids, boronic acids and acidic phenols including HOBT.

Cartridges format, Sorbent amount	#/box	SilactSPE SiSAX	SilactSPE Carbonate
1mL, 50mg	100	SiSAX-100.S.1.50	CO3-100.S.1.50
1mL, 100mg	100	SiSAX-100.S.1.100	CO3-100.S.1.100
3mL, 200mg	50	SiSAX-50.S.3.200	CO3-50.S.3.200
3mL, 500mg	50	SiSAX-50.S.3.500	CO3-50.S.3.500
6mL, 500mg	50	SiSAX-50.S.6.500	CO3-50.S.6.500
6mL, 1g	50	SiSAX-50.S.6.1000	CO3-50.S.6.1000
10mL LRC, 500mg	50	SiSAX-50.LRC.10.500	CO3-50.LRC.10.500
12mL, 2g	20	SiSAX-20.S.12.2000	CO3-20.S.12.2000
Reversible 0.7mL, 260mg	25	SiSAX-25.REV.1.260	CO3-25.REV.1.260
Reversible 2mL, 1000mg	25	SiSAX-25.REV.2.1000	CO3-25.REV.2.1000
96 well plate, 50mg	1	SiSAX-1.96W.50	CO3--1.96W.50
96 well plate, 100mg	1	SiSAX-1.96W.100	CO3-1.96W.100

**SilactSPE Florisil PR** (MgO<sub>3</sub>Si) : **Polar sorbent**

They present a basic character used to extract non-polar to moderately polar compounds from non-polar solvents.

They particularly suit for the retention of chlorinated pesticides, polychlorinated biphenyl (PCB's) and polysaccharides due to the magnesium ion.

**Product Information**

**Florisil PR type** : 150-200 μm

Cartridges format, Sorbent amount, #/box	SilactSPE Florisil PR
1mL, 50mg, 100/pk	FloPR-100.S.1.50
1mL, 100mg, 100/pk	FloPR-100.S.1.100
3mL, 200mg, 50/pk	FloPR-50.S.3.200
3mL, 500mg, 50/pk	FloPR-50.S.3.500
6mL, 500mg, 50/pk	FloPR-50.S.6.500
6mL, 1g, 50/pk	FloPR-50.S.6.1000
10mL LRC, 500mg, 50/pk	FloPR-50.LRC.10.500
12mL, 2g, 20/pk	FloPR-20.S.12.2000
Reversible 0.7mL, 300mg, 25/pk	FloPR-25.REV.1.300
Reversible 2mL, 900mg, 25/pk	FloPR-25.REV.2.900
96 well plate, 50mg, 1 unit	FloPR--1.96W.50
96 well plate, 100mg, 1 unit	FloPR-1.96W.100

**SilactSPE Dry**: contains sodium sulfate anhydrous (Na<sub>2</sub>SO<sub>4</sub>) in reversible cartridges.

Cart. format, Sorbent amount, #/box	SilactSPE Dry
Reversible 0.7mL, 800mg, 50/pk	Na2SO4-50.REV.1.800
Reversible 2mL, 2500mg, 50/pk	Na2SO4-50.REV.2.2500

**SilactSPE Na<sub>2</sub>SO<sub>4</sub>/Florisil** contains an upper layer of sodium sulfate anhydrous (Na<sub>2</sub>SO<sub>4</sub>) to dry the solution and a bottom layer of Florisil for the determination of hydrocarbons in water according to DIN-H-53/ ISO 9377-4.

Cart. format, Sorbent amount #/box	SilactSPE Na <sub>2</sub> SO <sub>4</sub> / Florisil
6mL PP, 2g+2g, 50/pk	FloNa2SO4-50.S.6.2g.2g
12mL PP, 3g+3g, 25/pk	FloNa2SO4-25.S.12.3g.3g

## SilactSPE Alumina-Acidic, Neutral and Basic

Alumina can present either cationic, neutral and acidic character. It is used in a similar fashion as for the SilactSPE Silica. The difference is that Alumina is more stable at high pH than silica.

**SilactSPE Alumina** particularly suit for the retention of aromatic compounds, aliphatic amines and compounds containing electronegative functions.

### Product Information

Alumina type : 60 Å, 0.9 g/mL, 50-200 µm

Cartridges format, Sorbent amount	#/box	SilactSPE Alumina Acidic	SilactSPE Alumina Neutral	SilactSPE Alumina Basic
1mL, 50mg	100	AluA-100.S.1.50	AluN-100.S.1.50	AluB-100.S.1.50
1mL, 100mg	100	AluA-100.S.1.100	AluN-100.S.1.100	AluB-100.S.1.100
3mL, 200mg	50	AluA-50.S.3.200	AluN-50.S.3.200	AluB-50.S.3.200
3mL, 500mg	50	AluA-50.S.3.500	AluN- 50.S.3.500	AluB- 50.S.3.500
6mL, 500mg	50	AluA-50.S.6.500	AluN-50.S.6.500	AluB-50.S.6.500
6mL, 1g	50	AluA-50.S.6.1000	AluN-50.S.6.1000	AluB-50.S.6.1000
10mL LRC, 500mg	50	AluA- 50.LRC.10.500	AluN- 50.LRC.10.500	AluB-50.LRC.10.500
12mL, 2g	20	AluA-20.S.12.2000	AluN-20.S.12.2000	AluB-20.S.12.2000
Reversible 0.7mL, 700mg	25	AluA-25.REV.1.700	AluN-25.REV.1.700	AluB-25.REV.1.700
Reversible 2mL, 2g	25	AluA- 25.REV.2.2000	AluN- 25.REV.2.2000	AluB- 25.REV.2.2000
96 well plate, 50mg	1	AluA-1.96W.50	AluN-1.96W.50	AluB--1.96W.50
96 well plate, 100mg	1	AluA-1.96W.100	AluN-1.96W.100	AluB-1.96W.100

## AttractSPE® Carbon

For the extraction of herbicides (EPA method 535)

A Graphitized Carbon Black sorbent. for absorption of pigments in food and small organic residues in water.

## AttractSPE® Carbon/Amine

For the cleanup of pesticides in food matrices prior to GC analysis

A two layer sorbents with Graphitized Black Carbon (GCB) and Aminopropyl modified silica sorbents

## AttractSPE® Carbon/PSA

For the cleanup of pesticides in food matrices prior to GC analysis

A two layer sorbents with Graphitized Black Carbon (GCB) and PSA modified silica sorbents

Product	Vol	Sorbent	25 cartridges/box	50 cartridges/box
<b>AttractSPE® Carbon</b>	6mL	500mg	Carb-25.S.6.500	Carb-50.S.6.500
<b>AttractSPE® Carbon/PSA</b>	3mL	250mg/250 mg	CarbPSA-25.S.3.250.250	CarbPSA-50.S.3.250.250
	6mL	500mg/500 mg	CarbPSA-25.S.6.500.500	CarbPSA-50.S.6.500.500
<b>AttractSPE® Carbon/Amine</b>	6mL	500mg/500 mg	CarbNH2-25.S.6.500.500	CarbNH2-50.S.6.500.500

# AttractSPE® on-line SPE



# ON-LINE SPE – DESCRIPTION AND PRODUCT LIST


## ON-LINE SPE PROCEDURE STEPS

SPE technique can be coupled on-line to HPLC for a high sensitivity or for limited amount of sample. An on-line SPE column containing the SPE sorbents is coupled. A three steps process is used:

**1- Sample injection:** The valve is configured with the injector directly in contact with the on-line SPE column. The sample is injected and goes through the on-line SPE column where the analytes remain.

**2- Washing:** a solution is used to wash out most interferences.

**3- Analysis:** The valve is switched. The analytes are eluted out of the sorbent by the LC mobile phase and transferred into the analytical column for their analyses.



Sample injection and washing (1 and 2)

Sample elution and analysis (step 3)

On-line SPE columns



Product	Product reference	Number column	I.D. (mm)	Lenght (mm)
<b>On-line AttractSPE® HLB columns</b>	OnlineSPE-HLB-1.2.20	1	2.1	20
	OnlineSPE-HLB-1.5.20	1	4.6	20
<b>On-line AFFINIMIP® PHENOLICS columns</b>	OnlineSPE-PHE-1.2.20	1	2.1	20
	OnlineSPE-PHE-1.5.20	1	4.6	20
<b>On-line AFFINIMIP® ESTROGENS columns</b>	OnlineSPE-EST-1.2.20	1	2.1	20
	OnlineSPE-EST-1.5.20	1	4.6	20

For other on-line SPE products, please contact us

AFFINISEP can provide you with on-line SPE of all products on demand.

# AttractSPE® SLE

Supported Liquid Extraction

# AttractSPE® LipRem

Proteins & lipid Removal

# AttractFiltration

Filtration

## SPE FOR REMOVAL OF PROTEINS & LIPIDS

### AttractSPE® LipRem

**For the removal of phospholipids of plasma sample**

**AttractSPE® LipRem** is a sorbent used for the removal of phosphorylcholine lipids from the plasma.

Cartridges format, Sorbent amount	#/box	AttractSPE® LipRem
1mL, 20mg	100	LipRem-100.S.1.20
3mL, 60mg	50	LipRem-50.S.3.50
6mL, 100mg	50	LipRem-50.S.6.100
96 well Plate	1	LipRem-1.96W.20
Reversible, 0.7mL, 100mg	25	LipRem-25.REV.1.F
	50	LipRem-50.REV.1.F

### SilactSPE Double fritted & SilactSPE Single fritted

**For the removal of proteins after precipitation**


**SilactSPE Double fritted & SilactSPE Single fritted** are cartridges with respectively one or two 20µm PE frits.

Cartridge volume	SilactSPE Double fritted 100 cartridges	SilactSPE Single fritted 100 cartridges
1mL	0-100.S.1.2F	0-100.S.1.1F
3mL	0-100.S.3.2F	0-100.S.3.1F
6mL	0-100.S.6.2F	0-100.S.6.1F
15mL	0-100.S.15.2F	0-100.S.15.1F
25mL	0-100.S.25.2F	0-100.S.25.1F
60mL	0-100.S.60.2F	0-100.S.60.1F
96 well plate – 1 unit	0-1.96W.2F	0-1.96W.1F

# AttractSPE®SLE FOR SUPPORTED LIQUID EXTRACTION

Supported Liquid Extraction (*a.k.a* SLE) is an alternative to LLE to pass from an aqueous media to an organic media without emulsion formation. This method is also useful to remove phospholipids and proteins from biological fluids.

AttractSPE®SLE contains an inert sorbent which absorbs water and enables the extraction of analytes with an organic solvent not miscible with water. This product advantageously replaces the phase transfer using liquid - liquid extraction and inherent problems such as emulsion formation. This process is easy to automatize, with a limited labour, glassware and organic solvent.



AttractSPE®SLE		
Loading volume	Format	50 cartridges/box
200µL	Column	SLE-50.S.200u
1mL	Column	SLE-50.S.1
2mL	Column	SLE-50.S.2
200µL	96 wellplate – 1unit	SLE-1.96W.200

# AttractFiltration FOR MEMBRANE FILTER CARTRIDGE

**AttractFiltration** is a filtration cartridges based on the use of a membrane to filtrate and remove particles with vacuum manifold (ACC-MAN1) or SPE automates before LC analysis.

A broad range of membranes is available and can enable a broad range of sample filtration. Available formats are 3mL, 6mL and 96 microfilter plate.

**AttractFiltration PES** with a PES membrane (hydrophilic, low protein binding) for water filtration

**AttractFiltration PTFE** with a PTFE membrane (hydrophobic, wide chemical compatibility, T resistance) for the filtraton of aggressive solutions

**AttractFiltration PVDF** with a PVDF membrane (hydrophobic, wide chemical compatibility, T resistance) for the filtraton of aggressive solutions

**AttractFiltration RC** with a Regenerated cellulose membrane (hydrophilic, solvent resistant, low non specific adsorption) for particle removal in solvents

**AttractFiltration Nylon** with a Nylon membrane (hydrophilic, high protein, RNA & DNA binding, high surface area) for a wide range of biological preparations

**AttractFiltration CA** with a Cellulose Acetate membrane (hydrophilic, low protein binding ) for protein filtration

Designation	Membrane	Pore size $\mu\text{m}$	3mL (100/box)	6mL (100/box)	96 filter plate – 1unit
<b>AttractFiltration PES</b>	PES	0.2	PES-100.S.3.2	PES-100.S.6.2	PES-1.96W.2
		0.45	PES-100.S.3.45	PES-100.S.6.45	PES-1.96W.45
<b>AttractFiltration PTFE</b>	PTFE	0.2	PTFE-100.S.3.2	PTFE-100.S.6.2	PTFE-1.96W.2
		0.45	PTFE-100.S.3.45	PTFE-100.S.6.45	PTFE-1.96W.45
<b>AttractFiltration PVDF</b>	PVDF	0.2	PVDF-100.S.3.2	PVDF-100.S.6.2	PVDF-1.96W.2
		0.45	PVDF-100.S.3.45	PVDF-100.S.6.45	PVDF-1.96W.45
<b>AttractFiltration RC</b>	Regenerated cellulose	0.2	RC-100.S.3.2	RC-100.S.6.2	RC-1.96W.2
		0.45	RC-100.S.3.45	RC-100.S.6.45	RC-1.96W.45
<b>AttractFiltration Nylon</b>	Nylon	0.2	NY-100.S.3.2	NY-100.S.6.2	NY-1.96W.2
		0.45	NY-100.S.3.45	NY-100.S.6.45	NY-1.96W.45
<b>AttractFiltration CA</b>	Cellulose acetate	0.2	CA-100.S.3.2	CA-100.S.6.2	CA-1.96W.2
		0.45	CA-100.S.3.45	CA-100.S.6.45	CA-1.96W.45

# Qcleanup


## Dispersive SPE QuEChERS and extraction salts



## Qcleanup - DESCRIPTION

QuEChERS is a sample pretreatment initially developed by Michelangelo Anastassiades and Steven Lehotay is mainly used for the analysis of multiple pesticides into fruits and vegetables. It is the acronym of **Quick, Easy, Cheap, Effective, Rugged and Safe**. Three main methods are currently used the original method, the European standardized method EN 15662 and the AOAC official method 2007.01. Each method required an extraction process with salts and a clean-up process with SPE cartridges or by dispersive SPE.

Schematic protocol of sample preparation with QuEChERS method



AFFINISEP supplies all products required to carry out QuEChERS according to AOAC or CEN including dispersive SPE products or SPE cartridges.

## Qcleanup DISPERSIVE SPE PRODUCTS

**Qcleanup** products for dispersive SPE are mixtures of powder in 2mL or 15mL centrifugation tubes for main scenarios encountered during pesticide analyses. This mixture contains magnesium sulfate anhydrous (MgSO<sub>4</sub>), primary secondary amine (PSA), carbon black (CB) or C18.

Method	Description	Nber/box	Product reference
<b>For General Fruits &amp; Vegetables</b>			
EN 15662	150mg MgSO <sub>4</sub> + 25mg PSA	100 tubes of 2mL	dSPE.EN.GFV.100.2
	900mg MgSO <sub>4</sub> + 150mg PSA	50 tubes of 15mL	dSPE.EN.GFV.50.15
AOAC 2007.01	150mg MgSO <sub>4</sub> + 50mg PSA	100 tubes of 2mL	dSPE.AOAC.GFV.100.2
	1200mg MgSO <sub>4</sub> + 400mg PSA	50 tubes of 15mL	dSPE.AOAC.GFV.50.15
<b>For Pigmented Fruits &amp; Vegetables</b>			
EN 15662	150mg MgSO <sub>4</sub> + 25mg PSA + 2.5mg CB	100 tubes of 2mL	dSPE.EN.PFV.100.2
	900mg MgSO <sub>4</sub> + 150mg PSA + 15mg CB	50 tubes of 15mL	dSPE.EN.PFV.50.15
AOAC 2007.01	150mg MgSO <sub>4</sub> + 50mg PSA + 50mg CB	100 tubes of 2mL	dSPE.AOAC.PFV.100.2
	1200mg MgSO <sub>4</sub> + 400mg PSA + 400mg CB	50 tubes of 15mL	dSPE.AOAC.PFV.50.15
<b>For Highly Pigmented and Fatty Fruits &amp; Vegetables</b>			
EN 15662	150mg MgSO <sub>4</sub> + 25mg PSA + 7.5mg CB	100 tubes of 2mL	dSPE.EN.HPFV.100.2
	900mg MgSO <sub>4</sub> + 150mg PSA + 45mg CB	50 tubes of 15mL	dSPE.EN.HPFV.50.15
AOAC 2007.01	150mg MgSO <sub>4</sub> + 50mg PSA + 50mg CB + 50mg C18	100 tubes of 2mL	dSPE.AOAC.HPFV.100.2
	1200mg MgSO <sub>4</sub> + 400mg PSA + 400mg CB + 400mg C18	50 tubes of 15mL	dSPE.AOAC.HPFV.50.15
<b>For Fatty and waxed Fruits &amp; Vegetables</b>			
EN 15662	150mg MgSO <sub>4</sub> + 25mg PSA + 25mg C18	100 tubes of 2mL	dSPE.EN.FWFV.100.2
	900mg MgSO <sub>4</sub> + 150mg PSA + 150mg C18	50 tubes of 15mL	dSPE.EN.FWFV.50.15
AOAC 2007.01	150mg MgSO <sub>4</sub> + 50mg PSA + 50mg C18	100 tubes of 2mL	dSPE.AOAC.FWFV.100.2
	1200mg MgSO <sub>4</sub> + 400mg PSA + 400mg C18	50 tubes of 15mL	dSPE.AOAC.FWFV.50.15

## Qcleanup EXTRACTION SALTS



**Qcleanup** extraction salts are the three main salts mixtures used in QuEChERS method.

QuEChERS methods	Description	Pouches / box	Product reference
Original method	4g $MgSO_4$ 1g NaCl	50	EXT.ORL.50
EN 15662	1g Trisodium citrate Dihydrate 0.5g Disodium hydrogencitrate sesquihydrate 1g NaCl and 4g $MgSO_4$	50	EXT.EN.50
AOAC 2007.01	1.5g Sodium Acetate and 6g $MgSO_4$	50	EXT.AOAC.50



2L - 500mL

- AttractSPE® Disks
- Passive samplers: POCIS, SPATT, Chemcatcher

## Environmental applications



# **AttractSPE<sup>®</sup> Disks for environmental applications**



**AttractSPE®Disks** are Solid Phase Extraction Disks for the extractions of a broad range of contaminants. **AttractSPE®Disks** are **thin, dense** and **uniform** SPE disks for retention of targeted analytes without any breakthrough. Our innovative SPE disks allow the best interactions with analytes and a maximal flow rates without any channeling.


### IMPORTANT FACT

#### Compared to other Disks available on the market :

- ✓ Higher density avoiding channeling
- ✓ No bubbling
- ✓ More diversity on chemistry **including HLB**
- ✓ More diversity in term of capacity (different particle sizes and porosities)
- ✓ Large diversity of applications for passive sampling (Chemcatcher and POCIS) as well as for protein and peptide purifications
- ✓ Can handle dirty samples
- ✓ Compatible with manual and automated holders
- ✓ Multimode or Specific application like disks for Glyphosate

AFFINISEP is manufacturer of its own resins from A to Z (from monomers to polymers) and can tailor-made chemistries with different physico-chemical properties.

- ❑ The membrane has a high exposed surface area of active particles
- ❑ AttractSPE® Disks are used to prepare samples of large volumes of water for environmental analysis (Compatible with many EPA methods)
- ❑ AttractSPE® Disks are available in broad diversity of Chemistry, HLB, C18, SDB, RPS, Anion and cation exchanges as well as multimode HLB and ionics
- ❑ Make possible the loading of high volume of water in a short time with excellent performances
- ❑ Usable as Passive sampler such as ChemCatcher
- ❑ Disks format for contaminants enrichment – 47mm or 90mm



## AttractSPE® Disks Environment

Designation	Description	Reference - 47mm -20/pk	Reference -90mm - 10/pk
<b>AttractSPE® Disks HLB</b>	HLB sorbent	SPE-Disks-HLB-47.T1.20	SPE-Disks-HLB-90.T1.10
<b>AttractSPE® Disks C18</b>	C18 sorbent	SPE-Disks-C18-47.T1.20	SPE-Disks-C18-90.T1.10
<b>AttractSPE® Disks C18 Polar</b>	C18 Polar sorbent	SPE-Disks-C18P-47.T1.20	SPE-Disks-C18P-90.T1.10
<b>AttractSPE® Disks SDB</b>	PS-DVB sorbent	SPE-Disks-DVB-47.T1.20	SPE-Disks-DVB-90.T1.10
<b>AttractSPE® Disks RPS</b>	Modified RPS sorbent	SPE-Disks-RPS-47.T1.20	SPE-Disks-RPS-90.T1.10
<b>AttractSPE® Disks Anion Exchange SR</b>	SAX sorbent	SPE-Disks-AN-47.T1.20	SPE-Disks-AN-90.T1.10
<b>AttractSPE® Disks Cation Exchange SR</b>	SCX sorbent	SPE-Disks-CAT-47.T1.20	SPE-Disks-CAT-90.T1.10
<b>AttractSPE® Disks Oil &amp; Grease</b>	For oil and grease	SPE-Disks-OIL-47.T1.20	SPE-Disks-OIL-90.T1.10
<b>AttractSPE® Disks Chelating</b>	For multivalent metal	SPE-Disks-MET-47.T1.20	SPE-Disks-MET-90.T1.10
<b>AttractSPE® Disks SAX-HLB</b>	SAX-HLB mixture	SPE-Disks-SAX-HLB-47.T1.20	SPE-Disks-SAX-HLB-90.T1.10
<b>AttractSPE® Disks High Spectrum</b>	Mix HLB-WCX-WAX for non-targeted screening	SPE-Disks-screening1-47.T1.20	SPE-Disks-screening1-90.T1.10
<b>AttractSPE® Disks High Spectrum</b>	Mix HLB-SCX-SAX for non-targeted screening	SPE-Disks-screening2-47.T1.20	SPE-Disks-screening2-90.T1.10
<b>AFFINIMIP® SPE Disks Picolinic Herbicides</b>	Based on <b>AFFINIMIP® SPE Picolinic herbicides</b> for extraction of picloram Clopyralid, Aminopyralid	SPE-Disks-PICO-47.T1.20	SPE-Disks-PICO-90.T1.10
<b>AFFINIMIP® SPE Disks Estrogens</b>	Based on <b>AFFINIMIP® SPE Estrogens</b>	SPE-Disks-EST-47.T1.20	SPE-Disks-EST-90.T1.10

## AttractSPE® disks Manifolds

One-, three- or six-station filtration manifolds allow the simultaneous extractions of several 1-L samples on a very simple and easy-to-handle way. The manifold is a very compact stainless steel device with a filtration glassware. Each station is controlled through an independent flow control valve.



Manifolds 47mm

Designation	Reference
AttractSPE® disks Manifolds - 1 station - 47mm	ACC-DISKSPE-G47-1
AttractSPE® disks Manifolds - 3 stations - 47mm	ACC-DISKSPE-G47-3
AttractSPE® disks Manifolds - 6 stations - 47mm	ACC-DISKSPE-G47-6
4L Polycarbonate Trap with 2x 1m vacuum hose	ACC-TRAP-4L



## AttractSPE® Prefilter Glassfiber

AttractSPE® Prefilter Glassfiber are required on top of AttractSPE® Disks to prevent clogging when loading water rich of suspended particles



Description	Reference	Description	Reference
25mm, 1µm, 50/pk	PF-GF-50.T1.25.1	90mm, 1µm, 50/pk	PF-GF-50.T1.90.1
25mm, 3µm, 50/pk	PF-GF-50.T1.25.3	90mm, 3µm, 50/pk	PF-GF-50.T1.90.3
47mm, 1µm, 50/pk	PF-GF-50.T1.47.1	100mm, 1µm, 50/pk	PF-GF-50.T1.100.1
47mm, 3µm, 50/pk	PF-GF-50.T1.47.3	100mm, 3µm, 50/pk	PF-GF-50.T1.100.3

# Passive Sampling Solutions

**POCIS**

**SPATT**

**Passive samplers using Disks  
Silicone rubber strips**

Passive sampling enables the monitoring of contaminants in water (surface water, groundwater, coastal water...) for a short (at least 7 days) to long period (with an average field deployment of one month) for which no power, maintenance and supervision is required. An average of the concentration of collected contaminants is measured in the laboratory.



## Advantages of Passive Sampling


- ❑ Can generate a time-weighted average (TWA) concentration of the contaminants in water
- ❑ Deployable in harsh conditions
- ❑ No a priori preparation or supervision
- ❑ Very simple use

### Passive Samplers and Applications

POCIS	SPATT	Chemcatcher & disks passive sampler	Silicon rubber
Polar Organic Chemical Integrative Sampler	Biotoxin tracking	Micropollutants and Metals	Hydrophobics compounds PAHs, PCBs, pesticides



### How is made it? Example of a POCIS



# Passive Sampling

## Steps from water to analysis

The sorbent collects the contaminant(s) in water. The organic compounds are then extracted from the sorbent in the POCIS, following a SPE procedure and analyzed using classical analytical methods by HPLC, and LC-MS, ....



## Broad range of applications for your sample preparation

**Hormons and EDCs**  
Natural & synthetic  
Estrogens  
Bisphenols & analogs  
Phenolics  
...

**Pesticides**  
Glyphosate & AMPA  
Aminopyralid,  
Clopyralid, Picloram  
Atrazine & derivatives  
DiuroPesticidesn  
...

**Drug residues**  
Carbamazepine  
Sulfamethoxazole  
Diclofenac  
Propranolol  
Erythromycin  
Tetracycline  
...

**Other contaminants**  
Caffeine  
PFOS, PFOA  
Biocides  
PAHs  
Biotoxins

### POCIS families

- **AttractSPE® POCIS HLB** for pharmaceutical drug residues
- **AttractSPE® POCIS Pesticides** for pesticides
- **AFFINIMIP® POCIS GLYPHOSATE** for Glyphosate and AMPA
- **AFFINIMIP® POCIS EDC** for the retention of Phenolic endocrine disruptors such as natural/synthetic estrogens, Bisphenols...

**Round POCIS for surface monitoring or rectangular for groundwater.**

#### Our POCIS with PRC

**Performance and reference compound (PRC)** is a compound not present in the environment (e.g., a deuterated molecule), which is spiked in the sorbent phase of the POCIS before its exposure as internal standards.

### Silicone rubber strips (SR)

Silicone rubber is a very powerful tool for the monitoring of non-polar contaminants such as PAHs, PCBs or pesticides.

Thanks to a very high exchange area, low limit of quantification of these contaminants can be reached.



#### Do you know ?

*AFFINISEP provides a complete range of SPE Disks suitable for **Chemcatcher and disks passive sampler**. use (see page SPE Disks for more information)*

### AttractSPE® SPATT Biotoxins

**AttractSPE® SPATT Biotoxins** is used to evaluate the contamination of shellfish with biotoxins in seawater such as Pectenotoxin, Yessotoxin, Okadaic acid/ Dinophysistoxin and Azaspiracids.

## Passive sampling PRODUCT LIST

Several kits are available **with or without a performance reference compounds (PRC)** to correct for in situ exposure known to affect uptake rates. These kits also include empty fritted cartridges to make easiest the extraction step of the contaminants.



Designation	Description	Composition	Reference
<b>AFFINIMIP® POCIS GLYPHOSATE</b>	POCIS containing AFFINIMIP® GLYPHOSATE - AMPA for the retention of glyphosate and AMPA	1 POCIS	POCIS.GLY.90.55.A .1
		Kit of 10 POCIS + empty fritted cartridges	POCIS.GLY.90.55.kit.10
<b>AFFINIMIP® POCIS EDC</b>	POCIS containing AFFINIMIP® Estrogens and AFFINIMIP® Bisphenols for the retention of endocrine disrupters such as natural/synthetic estrogens, Bisphenols...	1 POCIS	POCIS.EDC.90.55.A .1
		Kit of 10 POCIS + empty fritted cartridges	POCIS.EDC.90.55.kit.10
<b>AttractSPE® POCIS Pesticides</b>	POCIS containing mixture of sorbent for the retention of several pesticides	1 POCIS	POCIS.PEST.90.55.A.1
		Kit of 10 POCIS + empty fritted cartridges	POCIS.PEST.90.55.kit.10
		Kit of 1 POCIS with a sorbent containing DIA as PRC - 3 cartridges containing the sorbents with Désisopropylatrazine (DIA) d5 - empty fritted cartridges	POCIS.PEST.90.55.kit.1.DIA
<b>AttractSPE® POCIS HLB</b>	POCIS containing Attract HLB for the retention of pharmaceutical drug residues	1 POCIS	POCIS.HLB.90.55.A .1
		Kit of 10 POCIS + empty fritted cartridges	POCIS.HLB.90.55.kit.10
		Kit of 1 POCIS with a sorbent containing DIA as PRC - 3 cartridges containing the sorbents with Désisopropylatrazine (DIA) d5 - empty fritted cartridges	POCIS.HLB.90.55.kit.1.DIA

# Passive sampling PRODUCT LIST

## POCIS for groundwater

Designation	Definition	Reference
<b>AttractSPE® POCIS</b> <b>HLB for</b> <b>Groundwater</b>	1 POCIS containing Attract HLB for the retention of pharmaceutical drug residues & pesticides - <b>30cm x 5cm</b> - mass sorbent equivalent to 2 rounds POCIS	POCIS.HLB.30.5.A.1
	POCIS containing Attract HLB for the retention of pharmaceutical drug residues & pesticides - <b>1 unit = 3 pocis of 10cmx4cm equivalent to 2 round POCIS</b>	POCIS.HLB.10.4.A.3
<b>AFFINIMIP®SPE</b> <b>POCIS GLYPHOSATE</b> <b>for Groundwater</b>	1 POCIS containing AFFINIMIP® Glyphosate - <b>30cm x 5cm</b> - mass sorbent equivalent to 2 rounds POCIS	POCIS.GLY.30.5.A.1
	POCIS containing AFFINIMIP® Glyphosate in Groundwater - <b>1 unit = 3 pocis of 10cmx4cm equivalent to 2 round POCIS</b>	POCIS.GLY.10.4.A.3



### POCIS FOR GROUNDWATER

1 unit = **3 pocis of 10cmx4cm**



### POCIS FOR GROUNDWATER

1 unit = **30cm x 5cm**

AFFINISEP provides POCIS with specific shapes for groundwater.

Do not hesitate to contact for the POCIS of interest.

## Disks passive samplers

Designation	Description	Reference
<b>AttractSPE®Disks</b> <b>Passive Sampler HLB</b>	Disks - based passive samplers with AttractSPE®Disks HLB + ready to use PES membranes, 10/pk Outer diameter of 90mm and with a single-sided opening of 40mm. Compliant with 3-PS holders	DBPS.HLB.90.40.kit. 10
<b>AttractSPE® Disks</b> <b>Passive Sampler RPS</b>	Disks - based passive samplers with AttractSPE®Disks RPS + ready to use PES membranes, 10/pk Outer diameter of 90mm and with a single-sided opening of 40mm. Compliant with 3-PS holders	DBPS.RPS.90.40.kit. 10

Front and back sides of disks passive samplers- compliant with POCIS holder



## Passive sampling PRODUCT LIST

AttractSPE® SPATT			
Designation	Description	Composition	Reference
<b>AttractSPE® SPATT BIOTOXINS</b>	SPATT containing HP-20 sorbent for the retention of biotoxins. Nylon mesh membrane	10 SPATT Biotoxins	SPATT.BIOTOX.90.55 .300.A.10

Attract Silicone rubber		
Designation	Description	Reference
<b>Attract Silicone Rubber SR</b>	Silicone rubber strips – 10/pk – No PRC	SR.0.100.3.A1.10
	Silicone rubber strips with PRCs – 10/pk	SR.PRC.100.3.A1.10



AFFINISEP provides a complete set of accessories for the use of **Passive sampler**.

Designation	Description	Reference
<b>CANISTER – 3 POCIS</b>	1 Canister for 3 POCIS – Empty (holder not supplied)	CAN-3P.A.1
<b>CANISTER AND HOLDER FOR 3 POCIS</b>	1 Canister and 1 holder for 3 POCIS/DBPS	CH-3P.A.1
<b>CANISTER 24cm</b>	1 Canister 24cm – Empty - Requires two holders for 3 passive samplers each (not included)	CAN24.A.1
<b>CANISTER 29cm</b>	1 Canister 29cm – Empty - Requires 2 holders for 3 passive samplers each + 1 spider holder for silicone rubber/SPMD (not included)	CAN29.A.1
<b>HOLDER – 3 POCIS</b>	1 Holder for 3 POCIS	HOLD-3P.A.1
<b>HOLDER - SPIDER</b>	1 Holder – Spider for silicone rubber or SPMD	HOLD-SPI.A.1



**CANISTER – 3 POCIS**



**HOLDER – 3 POCIS**



**HOLDER – SPIDER**

**CANISTER - 29cm**  
Order on  
[www.affinisep.com](http://www.affinisep.com)

# **SPE ACCESSORIES**

**Vacuum Manifold**

**Mini Vap**

**Pump**

## AttractSPE<sup>®</sup> Vacuum Manifold

very flexible, allows you to control the flow and to process up to 12 or 24 samples simultaneously, to gain significantly time during sample preparation steps.



## Vacuum Manifold

**ACC-MAN2** Like all chromatography techniques, Use of SPE cartridges needs a precise control of flow rate for maintaining reproducible extractions. Solid Phase extraction Vacuum Manifold allows you to control the flow and to process up to 12 (12-port version) or 24 (24-port version) AFFINIMIP<sup>®</sup> SPE samples simultaneously, to gain significantly time during sample preparation steps.

## Mini PUMP

**ACC-PUMP** Diaphragm vacuum pump for solid phase extraction experiments

- 5.5L/min
- ~120 torr vacuum
- Oil-free
- Portable

## Vacuum pump trap

**ACC-TRAP-1L** SPE Vacuum pump trap kit

Installed between the manifold and the vacuum pump, it collects all liquids that are aspirated preventing contamination of the vacuum pump with a capacity of 1L.

## SPE ACCESSORIES

### SPE Adapter & Reservoir kit

**ACC-AR1** Tube adapters serve to pile one SPE tube on top of another to provide different selectivities. A larger empty syringe barrel can be stacked on top of a smaller SPE tube to act as a larger load reservoir. Or, they can serve as an adapter for positive pressure methods (e.g. from a syringe or air/ N2 line).

### Mini-Vap

**ACC-VAP1** The 6-Port Mini-Vap concentrator/evaporator processes six vials at one time. The Mini-Vap includes a needle valve for fine metering of air or nitrogen drying gas.

## SPE ACCESSORIES – Product list

SPE Accessories	Designation	Definition	Reference
Manifold	SPE Vacuum Manifold	12-port model	ACC-MAN2
SPE Adapter & Reservoir kit	SPE Adapter & Reservoir kit	Kit of 12 reservoirs 60ml and adapters for use with 1,3 & 6 mL cartridges	ACC-AR1
Mini-Vap	Mini Evaporator / Concentrator	6 port Mini-Vap Evaporator/Concentrator for use with 1 to 250mL containers	ACC-VAP1
Mini PUMP	Mini vacuum pump	Diaphragm vacuum mini pump, 5.5L/min	ACC-PUMP
Vacuum pump trap	SPE Vacuum pump trap kit	1L trap kit	ACC-TRAP-1L

# Examples of SPE applications



## EXAMPLES OF SPE APPLICATIONS

	ANALYTES	SPE product	MATRICES
	Multimycotoxins		
	Aflatoxins, Ochratoxin A, HT-2, T-2, Fumonisin, Zearalenone, Deoxynivalenol	<b>AFFINIMIP® SPE</b> <b>Multimyco LCMSMS</b>	Cereals
	Fumonisin AND Zearalenone	<b>AFFINIMIP® SPE</b> <b>FumoZON</b>	Maize, Maize-based baby food
	Single Mycotoxin		
	Patulin	<b>AFFINIMIP® SPE Patulin</b>	All Apple-based products (Juice, puree, concentrate...)
	Zearalenone	<b>AFFINIMIP® SPE</b> <b>Zearalenone</b>	Maize, Wheat, Cereal-based baby food, Rice, Edible corn oil
	Ochratoxin A	<b>AFFINIMIP® SPE</b> <b>Ochratoxin A</b>	Wheat, Maize, Pepper, Paprika, Red and White Wine
	Deoxynivalenol (DON)	<b>AFFINIMIP® SPE</b> <b>Deoxynivalenol</b>	Wheat, Maize, Oat
<b>Endocrine Disruptor</b>	Estrone, 17 $\alpha$ -Estradiol, 17 $\beta$ -Estradiol, Estriol, 17 $\alpha$ -Ethinylestradiol	<b>AFFINIMIP® SPE</b> <b>Estrogens</b>	Water, Serum, Plasma
	Bisphenol A, Bisphenol AP, Bisphenol AF, Bisphenol B, Bisphenol S, Bisphenol F...	<b>AFFINIMIP® SPE</b> <b>Bisphenols</b>	A broad variety of liquid and solid foods
	Parabens	<b>AFFINIMIP® SPE</b> <b>Phenolics</b>	Shampoo, cream
	Phenolic compounds	<b>AFFINIMIP® SPE</b> <b>Phenolics</b>	Food matrices
	Bisphenols & Alkyl phenols	<b>AttractSPE®Disks C18</b>	Water
	Endocrine disruptors	<b>AttractSPE®Disks HLB</b>	Water
<b>Drug Residues</b>	Amphetamine, Methamphetamine, MDA, MDMA, MDEA	<b>AFFINIMIP® SPE</b> <b>Amphetamines</b>	Serum, Urine
	Zeranol, Zearalanone, $\alpha$ and $\beta$ Zearalanol, $\alpha$ and $\beta$ Zearalenol, Resorcylic acid lactones	<b>AFFINIMIP® SPE</b> <b>Zeranol Residues</b>	Urine, Plasma
	Chloramphenicol	<b>AFFINIMIP® SPE</b> <b>Chloramphenicol</b>	Honey, Urine, Shrimp
	Tamoxifen	<b>AFFINIMIP® SPE</b> <b>Tamoxifen</b>	Urine

See our application notebook for more applications and details.

## EXAMPLES OF SPE APPLICATIONS

	ANALYTES	SPE product	MATRICES
Antibiotics and Drugs residues	Nicotine, Procainamide	<b>AttractSPE® HLB</b>	Urine
	Caffeine	<b>AttractSPE® HLB</b>	Urine, Water
	Propranolol	<b>AttractSPE® HLB</b>	Urine, Water
	<b>Tetracyclines</b> - Tetracycline, Oxytetracycline, Chlortetracycline, Doxycycline	<b>AFFINIMIP® SPE Tetracyclines</b>	Milk
	<b>Sulfonamides</b> – Sulfadimethoxine , Sulfamethoxypridazine...	<b>AttractSPE® SCX</b>	Milk
	Caffeine, Acetaminophen, Diclofenac, Ibuprofen, Ketoprofen, Naproxen, Carbamazepine	<b>AttractSPE® HLB</b>	Waste water, water
	<b>Antibacterial Aminoglycosides</b> - Streptomycin, Dihydrostreptomycin,...	<b>AttractSPE® HLB</b>	Tissue, Milk
	<b>Antibiotics</b> – Quinolones, Macrolides, Lincosamides, Sulfonamides, Penicillins, Cephalosporine, Pleuromutilins, Diamino pyrimidine derivatives	<b>AttractSPE® HLB</b>	Tissue, Milk
	<b>NSAID (Non Steroidal Anti inflammatory drug)</b> - Salicylic acid, Phenylbutazone, Flunixin, Tolfenamic acid, Meloxicam, Desoximetasone, Ketoprofen	<b>AttractSPE® HLB</b>	Tissue
	<b>Penicillin based antibacterials</b> - Ampicillin, Amoxicillin...	<b>AttractSPE® HLB</b>	Tissue
	<b>Glucocorticoids</b> - Cortisone, Corticosterone, Aldosterone, Betamethasone, Dexamethasone, Flumethasone, Prednisone, Prednisolone, Methylprednisolone	<b>AttractSPE® HLB</b>	Tissue
	Erythromycin and Clindamycin	<b>AttractSPE® HLB</b>	Tissue
	Praziquantel and Tiamulin	<b>AttractSPE® HLB</b>	Tissue
	Cephalexin	<b>AttractSPE® HLB</b>	Fish
	Quinoxaline-2 -carboxylic acid and 3-methyl quinoxaline-2-carboxylic acid	<b>AttractSPE® SAX</b>	Muscle, Liver, Kidneys
	Vancomycin	<b>AttractSPE® SCX</b>	Fish
	Valnemulin and Tiamulin	<b>AttractSPE® HLB</b>	Fish
	Phenolic compounds	<b>AFFINIMIP® SPE Phenolics</b>	biological matrices
	Hormons, sex steroids, Pharmaceutical compounds	<b>AttractSPE® Disks HLB</b>	Water
	Pharmaceutical compounds and analytes containing carboxylic acid groups	<b>AttractSPE® Disks Anion Exchange - SR</b>	Water

## EXAMPLES OF SPE APPLICATIONS

	ANALYTES	SPE product	MATRICES
Pesticides – Herbicides - Biocides	Glyphosate, AMPA	<b>AFFINIMIP® SPE Glyphosate</b>	Water
	Aminopyralid, Clopyralid, Picloram	<b>AFFINIMIP® SPE Picolinic Herbicides</b>	Water, Compost, Soil
	<b>16 common pesticides</b> - Linuron, Iprodione, Desisopropylatrazine, Desethylatrazine, Aldocarb, Simazine, Carbofuran, Metalaxyl, Atrazin, 2, 4-D, Metazachlor, Dicloran, Phenmedipham, Procymidone, Fenitrothion, Vinclozolin	<b>AttractSPE® HLB</b>	Water
	<b>Triazine Herbicides</b> - Simazine, Cyanazine, Atrazine...	<b>AttractSPE® HLB</b>	Water
	<b>Acetamide Herbicides</b> - Metolachlor and metabolites, Alachlor...	<b>AttractSPE® HLB</b>	Water
	<b>Fungicides</b> - Carbendazim, Thiabendazole	<b>AttractSPE® SCX</b>	Fruit Juice
	<b>Organotin compounds</b>	<b>SilactSPE Organotin</b>	Water, soils
	<b>Pesticides by GC-MS</b> : Metamidophos, Dichlorvos, Acephate, Trifluralin, Diazinon, Chlorothalonil, Dimethipin, Vinclozoline, Methyl parathion, Methyl primophos, Triadimenol-1, DDE, Cypermethrin-3, Difenconazole-1, Imibenconazole, Tebuthiuron, Bromacil...	<b>AttractSPE® Carbon/PSA</b>	Food matrices
	Benzidine & Nitrogen-containing Pesticides	<b>AttractSPE® Disks C18</b>	EPA 553.1
	Chlorinated acids	<b>AttractSPE® Disks XC</b>	EPA 515.2
	Endothall	<b>AttractSPE® Disks HLB AttractSPE® Disks C18 AttractSPE® Disks Anion Exchange - SR</b>	EPA 548; EPA 548.1 Rev. 1
	Nitrogen- & Phosphorus-Containing Pesticides	<b>AttractSPE® Disks C18</b>	EPA 507
	Organochlorine Pesticides	<b>AttractSPE® Disks HLB AttractSPE® Disks C18</b>	EPA 8081. EPA 608
Organophosphorus Pesticides	<b>AttractSPE® Disks C18</b>	EPA 1614; EPA 1657	

See our application notebook for more applications and details.

## EXAMPLES OF SPE APPLICATIONS

	ANALYTES	SPE product	MATRICES
<b>PAHs – POPs – Perfluorinated compounds</b>	<b>Perfluorinated compounds</b> (Perfluorobutanoic acid (PFBA), Perfluoropentanoic acid (PFPeA), Perfluorohexanoic acid (PFHxA), Perfluoroheptanoic acid (PFHpA), Perfluorooctanoic acid (PFOA), Perfluorononanoic acid (PFNA), Perfluorodecanoic acid (PFDA), Perfluorotetradecanoic acid (PFTA), Perfluorobutanesulfonic acid (PFBS), Perfluorohexane sulfonic acid (PFHxS), Perfluorooctanesulfonic acid (PFOS))	<b>AttractSPE® PFAS</b>	Water
	<b>Hydroxylated Polycyclic Aromatic Hydrocarbons</b> - 2-Naphtol, 2-Hydroxyfluorene, 9-Phenanthrol...	<b>AFFINIMIP® SPE Phenolics</b>	Contaminated soils
	PBDEs, Dioxins & Furans, PAHs	<b>AttractSPE® Disks HLB</b> <b>AttractSPE® Disks C18</b>	Water, EPA 1613; EPA 625
	PAHs	<b>AttractSPE® Disks C18</b>	EPA 550.1
	PCBs	<b>AttractSPE® Disks HLB</b> <b>AttractSPE® Disks C18</b>	EPA 8082; EPA 1668
	<b>Polycyclic Aromatic Hydrocarbons (PAH)</b>	<b>AFFINIMIP® SPE PAHs</b>	Fats and oil
<b>AttractSPE® HLB</b>		Waste water	
<b>SilactSPE CN/SiOH</b>		soil	
<b>Phenolics</b>	Guaiacol	<b>AFFINIMIP® SPE Phenolics</b>	Wines, water
	Carnosic acid	<b>AFFINIMIP® SPE Phenolics</b>	Meat, water
	Hydroquinone	<b>AFFINIMIP® SPE Phenolics</b>	Water

See our application notebook for more applications and details.

## EXAMPLES OF SPE APPLICATIONS

	ANALYTES	SPE product	MATRICES
<b>Ions removal</b>	<b>Transition metals ions</b>	<b>AttractSPE® IDA</b>	Aqueous solution
	Removal of anionic contaminants and neutralization of highly acidic samples	<b>AttractSPE® SAX-HCO<sub>3</sub></b>	Aqueous solutions
	Removal of alkaline earth and neutralization of basic samples	<b>AttractSPE® PS-H</b>	Aqueous solutions
	Removal of Halides ions (chloride, iodide, bromide)	<b>AttractSPE® PS-Ag</b>	Aqueous solutions
	Removal of sulfate ions	<b>AttractSPE® PS-Ba</b>	Aqueous solutions
	Removal of WATER – drying	<b>SilactSPE Dry</b>	Organic solutions
<b>Biological application</b>	Removal of phospholipids	<b>AttractSPE® LipRem</b>	plasma
	Removal of precipitated proteins	<b>SilactSPE Double fritted &amp; Single fritted</b>	Aqueous solutions
	Supported liquid extraction	<b>SilactSPE SLE</b>	Aqueous solutions
	NNAL (biomarker of smokers)	<b>AFFINIMIP® SPE NNAL</b>	Urine
	Dopamine, Noradrenaline, Adrenaline, ...	<b>AFFINIMIP® SPE Catecholamines</b>	Plasma, Serum
	Metanephrine, Normetanephrine and 3-Methoxytyramine, ...	<b>AFFINIMIP® SPE Metanephrines</b>	Plasma, Serum
	Fractionation and desalting of peptides and proteins	<b>AttractSPE® Disks Tips – Stagetips</b> <b>AttractSPE® Disks Spin</b> <b>AttractSPE® Disks 96 plate</b> <b>AttractSPE® Disks Cartridge</b>	Proteins digest, serum, urine, biological fluids
<b>Miscellaneous</b>	Melamine	<b>AttractSPE® SCX</b>	Milk, food
	Cyanuric acid	<b>AttractSPE® SAX</b>	Milk
	<b>ARTIFICIAL SWEETENERS</b> - Acesulfame, Aspartame, Cyclamate, Neohesperidine dihydrochalcone, Saccharin, Sucralose	<b>AttractSPE® HLB</b>	Water
	<b>COCAINE AND MAIN METABOLITES</b> - Cocaine, benzoylecgonine and ecgonine methyl ester	<b>AttractSPE® HLB</b>	Waste water
	<b>Unknown and not identified contaminants and metabolites</b>	<b>AttractSPE® Large Spectrum</b> (Non- target screening)	Waste water
	<b>Hydrocarbons</b> in water (ISO9377-4)	<b>SilactSPE Na<sub>2</sub>SO<sub>4</sub> / Florisil</b>	water

## EXAMPLES OF AttractSPE® Disks APPLICATIONS

	ANALYTES	EPA METHODS	SPE product
Extraction of contaminants in water with high flow rate	Acids and Base/Neutrals including PCBs	EPA 625	AttractSPE® Disks HLB
	Amines		AttractSPE® Disks Cation Exchange - SR
	Benzidine & Nitrogen-containing Pesticides	EPA 553.1	AttractSPE® Disks C18
	Carbonyl Compounds	EPA 8315	AttractSPE® Disks HLB AttractSPE® Disks C18
	Carbonyl Compounds & Formaldehyde	EPA 554	AttractSPE® Disks C18
	Chlorinated acids	EPA 515.2	AttractSPE® Disks SDB
	Chlorinated Pesticides	EPA 508.1	AttractSPE® Disks C18
	Dioxins & Furans	EPA 1613	AttractSPE® Disks C18
	Endothall	EPA 548; EPA 548.1 Rev. 1	AttractSPE® Disks HLB AttractSPE® Disks C18 AttractSPE® Disks Anion Exchange - SR
	Explosives Residues (HDX, RDX)		AttractSPE® Disks RPS
	Haloacetic Acids and Dalapon	EPA 552.1 Rev.1	AttractSPE® Disks Anion Exchange - SR
	Nitroaromatics & Nitramines	EPA 8330	AttractSPE® Disks HLB
	Nitrogen- & Phosphorus-Containing Pesticides	EPA 507	AttractSPE® Disks C18
	N-Methylcarbamate	EPA8318	AttractSPE® Disks HLB
	Oil & grease		AttractSPE® Disks Oil & Grease
	Organic Compounds	EPA 525	AttractSPE® Disks C18
	Organochlorine Pesticides	EPA 8081. EPA 608	AttractSPE® Disks HLB AttractSPE® Disks C18
	Organophosphorus Pesticides	EPA 1614; EPA 1657	AttractSPE® Disks C18
	PAHs	EPA 550.1	AttractSPE® Disks C18
	PCBs	EPA 8082; EPA 1668	AttractSPE® Disks HLB AttractSPE® Disks C18
	Phthalate Esters	EPA 8061	AttractSPE® Disks C18
	Phthalate & Adipate Esters	EPA 506	AttractSPE® Disks C18
	Phenylurea compounds	EPA 532	AttractSPE® Disks HLB AttractSPE® Disks C18
Semivolatile Organic Compounds	EPA 8270	AttractSPE® Disks HLB	
Hormons, sex steroids, PAHs, PPCPs, Pharmaceutical compounds, Endocrine disruptors		AttractSPE® Disks HLB	
Bisphenols & Alkyl phenols, PBDEs, Dioxins & Furans, Phthalates, Herbicides, PAHs, Carbaryl, Microcystins		AttractSPE® Disks C18	
Pesticides, Pharmaceutical compounds and analytes containing carboxylic acid groups		AttractSPE® Disks Anion Exchange - SR	

# CUSTOM-MADE PRODUCTS & PRIVATE LABEL SERVICES

- Development of custom products
- Private labelling of existing products
- Customization of existing products



Food / Feed Safety



Environment



Cosmetics




Life Science

**AFFINISEP offers full services for the design and synthesis of polymers complying with your specifications.**

With fully integrated technology platforms in polymer and analytical chemistry, AFFINISEP has been partner in more than 56 projects and helps its customers by innovative solutions for their complexes challenges.

AFFINISEP has developed a library of monomers giving a family of selective stationary phases based on its proprietary technology, which have shown a strong potential for the extraction and purification of various compounds.



### Why trusting AFFINISEP?

#### AFFINISEP supplies

- The most comprehensive bank of sorbents
  - From Silica to Polymers
  - From conventional to very selective
  - MIPs, Reversed-phase, ionic exchange, etc...
- Our Expertise on MIP, sample preparation, SPE protocol and detection kit development
- Quick and efficient development
- Reactivity

# PROJECT DEVELOPMENT

No matter if we run a short term project (2-3 days) or long term project ( 4 to 6 months) we always follow a well-established procedure. The following scheme describes an example of a procedure for the development of a custom-made product based on our customer's requirements.


After signing of Confidential Disclosure Agreements (CDA), key data are exchanged. A quotation based on our experience in separation science including a process sheet is developed. During the whole procedure, a permanent feedback is established with you.

## Procedure for custom-made polymer phase

**Step 1:** We offer you a screening of our library which consists of several hundreds AFFINISEP phases to find the suitable one for your separation problem. The knowledge of the structure of the target substance, its functional groups and the solubility data give us valuable hints for the choice of the screening phases.

**Step 2:** For the selected phase, a protocol is implemented for your application. Then samples and / or prototypes are delivered to you for testing, evaluation and validation.

**Step 3:** When the selected phase suits your application and has been validated, a scale-up is planned. A QC report is delivered with the product. The format of the product is correlated to your application and can be bulk material, SPE cartridges, HPLC columns etc...




If you need the development of new polymer for your application, please send us an email to [contact@affinisep.com](mailto:contact@affinisep.com) or describe your needs using <https://www.affinisep.com/technical-and-customer-support/technical-support/> You can describe your application and our scientists will shortly evaluate your queries before contacting you as soon as possible.

## PRIVATE LABEL SERVICES

**AFFINISEP offers a range of options that enable you to take advantage of our technology and production capacity to provide the latest products into your own portfolio.**

We work with you to provide tailored solutions that meet your specific requirements:






### Why trusting AFFINISEP?

- Significant R&D savings using our ready to use kits
- Benefit from our Expertise Quick and efficient development
- Dedicated technical support
- Ensure product quality made in France
- Ensure product availability by buying from manufacturer
- Ensure confidentiality

## THE ART OF MAKING SAMPLE PREPARATION EASIER

### About

AFFINISEP is a **worldwide expert in sample preparation applications**. AFFINISEP is dedicated to the development of analytical applications in various fields such as water, biological fluids, food, feed analysis and proteomics with a complete set of products and services for sample preparation.

 Brands	 Applications	 Matrices
AFFINIMIP SPE® AttractSPE® SilactSPE ...	Sample Preparation Passive Sampling Filtration Microextraction of peptides/proteins	Food, Feed, Soil, Oil, Water, Biological fluids, Proteolytic digestion

**Analytical chemists can find any solution for sample preparation, selective extraction and sample clean-up needs in various sectors: food and feed safety and quality, life science and quality control, clinical diagnosis, environment and doping.**

**In addition, proteomics users can find a complete set of microelution products for protein/peptides fractionation or desalting.**

### ORDERING INFORMATION



For any order, please, choose one of the following ways:

❖ On our website : [www.affinisep.com](http://www.affinisep.com)

Find all our references and benefit from a **FAST** and **SIMPLE** process !

❖ By email:

❖ [contact@affinisep.com](mailto:contact@affinisep.com)

#### Information about new products or new matrices:

❖ **Newsletters**

Stay up-to-date on AFFINISEP's latest innovations by subscribing to our newsletter, delivering new information straight to your inbox ! Register on our website.