




AFFINISEP

Every Micro-elution challenge has its
AttractSPE™ Disks Solution

AttractSPE™ Disks

96-WELL SPE PLATES






AttractSPE™ Disks Tips are great choice for peptide desalting & fractionation, proteomics, biomarker discoveries and biological applications



ADVANTAGES

- Capacity : 10-200 µg
- A 96 well tip shape for a high throughput screening
- Various chemistries C18, C4, C8, DVB, Ion exchange (WCX, SAX, WAX, SCX), HLB, RPS, and multilayers mixtures.
- Peptides and proteins purification, fractionation and desalting, proteins precipitation, drugs and metabolites screening



<https://www.affinisep.com>
contact@affinisep.com

